



**GORE® VIABAHN®**

Endoprosthesis  
with PROPATEN Bioactive Surface\*

OPEN MORE  
POSSIBILITIES



\* As used by Gore, PROPATEN Bioactive Surface refers to Gore's proprietary CBAS® Heparin Surface.

***Together, improving life***

1996

Original GORE® HEMOBAHN® Endoprosthesis introduced in Europe

2008

GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface approved in Europe

5–8 mm devices decreased in profile by one French size

2003

TIP to HUB deployment introduced on 6–8 mm devices


2009

Laser technology enables the new contoured edge at proximal end

9–13 mm devices introduced with 0.035" guidewire compatibility

2010

25 cm length:  
Longest stent-graft introduced in Europe



# Continued innovation for durable outcomes and unmatched versatility

The GORE® VIABAHN® Device is a leader among stent grafts. Decades of partnership with clinicians around the globe has resulted in unparalleled performance across multiple indications.\*

2011


GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface 5–8 mm devices decreased in profile by one French size

2016

Radiopaque markers introduced on 5–8 mm devices

2022

GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface 7.5 cm length introduced for the 5–9 mm devices



2014

Receives CE mark for the treatment of symptomatic venous obstruction

2021

GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface: up to 3 French size profile reduction and addition of radiopaque markers on 9–13 mm devices

\* GORE® VIABAHN® Endoprosthesis. W. L. Gore & Associates website. Accessed August 5, 2022. <https://www.goremedical.com/eu/VIABAHN/references>

# Backed by a growing body of clinical data in various clinical indications

The GORE® VIABAHN® Device has become a go-to device for physicians' most challenging cases



OVER  
**2,300**  
scientific publications†



OVER  
**5,650**  
hospitals using GORE® VIABAHN® Device†



OVER  
**1,000,000**  
implanted in patients worldwide\*†

Cephalic arch and central vein stenosis


Isolated visceral artery aneurysms

Dialysis access circuit stenosis

Iliac artery stenosis

(In-stent) restenosis of the superficial femoral and proximal popliteal artery

Popliteal artery aneurysms



\* GORE® VIABAHN® Endoprosthesis. W. L. Gore & Associates website. Accessed August 5, 2022. <https://www.goremedical.com/eu/VIABAHN/references>  
† Data shown is representative of all generations of the GORE® VIABAHN® Device.

# Iliac artery stenosis

The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface can be used for treating lesions in the iliac arteries.

High primary patency even in the most challenging disease: Demonstrated 87% primary patency at three years.<sup>1</sup>

Durable clinical study outcomes in complex cases:


Stent grafts have demonstrated superiority over bare metal stents (BMS) for treating complex iliac lesions.<sup>1,2</sup>

Self-expanding stent grafts, at three years, have demonstrated improved patency over BMS when treating TASC D iliac lesions.<sup>1</sup>



Before: Flush ostial CTO of the right common iliac artery with reconstitution of the proximal common femoral artery.

Higher patency out to three years for TASC II D lesions<sup>1</sup>:



Three-year self-expanding covered stents versus bare metal stents in iliac artery occlusions; primary patency rates at three years



After: Post-placement of 7 mm x 150 mm GORE® VIABAHN® Device and 7 mm x 59 mm balloon expandable covered stent.

Images courtesy of Barry Weinstock, M.D. Used with permission.

# Arteriovenous (AV) Access

Proven success in AV Access, even the most challenging cases, including:


- Early percutaneous transluminal angioplasty (PTA) failures
- Lesions at points of flexion
- Thrombosed grafts

## High primary patency even in the most challenging disease:


Increased trend in primary patency in thrombosed grafts of both the target lesion and the circuit by ~50% when compared to PTA at six months.<sup>3</sup>

Provided consistent patency independent of the number of times a lesion has previously been treated.<sup>3</sup>

## Durable treatment of thrombosed AV grafts



Before: Graft thrombosis secondary to outflow stenosis at the venous anastomosis of an AV graft.




After: At 60 months post-placement, the GORE VIABAHN Device has maintained secondary patency without any further episodes of thrombosis.

---

Durable clinical study outcomes in complex cases:  
83% access secondary patency and zero device  
fractures at two years when placed across the elbow.<sup>4</sup>

---



Proven to reduce reinterventions:  
Lowered mean number of  
interventions over two years  
in thrombosed grafts when  
compared to PTA<sup>5</sup>

### Recommendations for optimal outcomes in AV Access:

- Outflow wall apposition to the outflow vein is not necessary for quality outcomes
- Follow the IFU recommendation for 5–20% oversizing using the graft inner diameter as the target vessel<sup>§</sup>
- Do not use PTA outside of the device<sup>§</sup>

**A post-hoc analysis of the Gore REVISE Clinical Study demonstrated that undersizing the stent graft to the outflow vein trended toward increased patency<sup>\*,6</sup>**

6-month outcomes	Device apposition relative to the outflow vein	
	Undersized <sup>†</sup>	Apposed <sup>‡</sup>
Target lesion primary patency	60%	47%
Circuit primary patency	47%	40%

Note: The GORE® VIABAHN® Device should always be sized 5% to 20% greater than the AV graft diameter per the *Instructions for Use*.

Images courtesy of Daniel V. Patel, M.D. Used with permission.

\* Caution should be used when interpreting post-hoc analysis.

† The difference between the diameter of the vein and the device is  $\geq 1$  mm.

‡ The difference between the diameter of the vein and the device is  $< 1$  mm.

§ Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com)


# Superficial femoral artery (SFA)

Strong clinical performance in the most challenging cases including long SFA lesions and chronic total occlusions (CTOs).


## High primary patency even in the most challenging disease:

- 88% 12-month primary patency in SFA lesions averaging 22 cm in length<sup>7</sup>

Long SFA lesion of the right SFA



Before: Proximal SFA disease and mid-SFA occlusion.



After: Post-placement of three 5 mm GORE<sup>®</sup> VIABAHN<sup>®</sup> Devices.

Proven patency for complex SFA lesions across seven multicenter, prospective, randomized or single arm studies<sup>7-13</sup>

1,089

Limbs studied

71%

CTOs<sup>†</sup>

23 cm

Average lesion length<sup>†</sup>

80%

Average primary patency<sup>‡</sup>

Images courtesy of James Persky, M.D. Used with permission.

\* Data on file 2020; W. L. Gore & Associates, Inc; Flagstaff, AZ.


† Weighted average lesion length.

‡ One-year weighted average primary patency.




## Durable clinical study outcomes in complex cases:

Comparable clinical results to above the knee surgical bypass (both prosthetic<sup>14</sup> and native vein<sup>10</sup>).



Three-year secondary patency  
in complex disease

(27 cm average lesion length, 93% CTOs)<sup>15</sup>



Five-year freedom from target  
lesion revascularization (fTLR)<sup>16</sup>

## Recommendations for optimal outcomes in the SFA

### Device sizing considerations:

- Treat all of the disease (stent “normal to normal”)<sup>15</sup>
- Overlap devices by at least 1 cm<sup>8</sup>
- Avoid excessive oversizing (> 20%)<sup>15</sup>

### Implantation considerations:

- Ensure adequate inflow and outflow<sup>15</sup>
- Post dilate<sup>8</sup>


- Do not use PTA outside of the device<sup>15</sup>
- Place device at the SFA origin if proximal SFA disease is present<sup>15</sup>

### Follow-up considerations:


- Regular duplex ultrasonography follow-up<sup>17</sup>
- Prescribe appropriate antiplatelet therapy<sup>15</sup>
- Treat progressing disease<sup>17</sup>

# In-stent restenosis (ISR) of the SFA

Durable treatment for complex in-stent restenotic lesions.



Before: Diffuse SFA ISR in long-stented segment in the SFA.



After: Completion angiogram after placement of GORE® VIABAHN® Devices for ISR in the SFA.

## High primary patency even in the most challenging disease:

75% one-year primary with an average lesion length of over 17 cm.<sup>18</sup>

Fewer than one third the number of patients required an intervention at one year compared to PTA.<sup>15</sup>


## Durable clinical study outcomes in complex cases:

Four times greater primary patency compared to PTA at two years.<sup>18</sup>

More than three times greater fTLR compared to PTA at two years.<sup>15</sup>


## Proven to reduce reinterventions:

Fewer patients had reintervention procedures compared to PTA at two years.<sup>15</sup>



primary patency at one year<sup>18</sup>

17.3 cm mean lesion length<sup>18</sup>



fTLR at one year<sup>18</sup>

#### Recommendations for optimal outcomes in ISR:

- Extend the device at least 1 cm proximally and distally from the previously placed stent<sup>15</sup>
- Cover any adjacent disease at least 1 cm beyond the proximal and distal margins of the lesion<sup>15</sup>
- Follow the IFU recommendation for 5–20% oversizing using the healthy vessel diameter immediately proximal and distal to the lesion<sup>15</sup>
- Ensure guidewire has traversed the lesion intraluminally before completing PTA<sup>15</sup>

# Popliteal Artery Aneurysm (PAA)

Clinical performance in challenging PAA cases

Endovascular repair of popliteal aneurysms is associated with acceptable long-term patency and a very low risk of limb loss.<sup>19</sup>

## Endovascular Repair vs. Open Repair

### Early Outcomes<sup>20</sup> (30 days)

Hospital stay ↓  
Wound complications ↓  
Primary Patency 1Y ↓

### 3Y Outcomes<sup>20</sup>

Primary patency ↔  
Secondary patency ↔  
Mortality ↔  
Amputation ↔

↓ Lower than open repair

↔ No statistical difference to open repair

## Endovascular Repair – Long-Term Outcomes

### 10Y Outcomes<sup>19</sup>

51% Primary patency  
60% Secondary patency  
79% Amputation free survival  
0% Major amputations




Figure 1.  
Large (5 cm diameter) left popliteal artery aneurysm



Figure 2a.  
GORE® VIABAHN® Device placement to exclude the aneurysm





Figure 2b.  
Completion angiography showing good conformability to the artery tortuosity and total exclusion of the aneurysm by GORE® VIABAHN® Device

Images courtesy of Barry S. Weinstock, MD. Used with permission.

# Technology and clinical benefits

The unique design of the GORE® VIABAHN® Device enables treatment of even the most challenging peripheral cases.



## Performs as an endoluminal bypass:

- Covers and excludes diseased and irregular tissue.
- Provides a barrier from tissue ingrowth, minimizing ISR.

## Conformable yet durable design:

- Like with all Gore single nitinol wire stents, the design and frame construction reduces strain to provide mechanical durability.
- Proven flexibility maintains flow at points of flexion and increases anatomical options.

## Ease of use:

- Robust configurations cover a broad range of patient needs.
- Radiopaque markers enhance endoprosthesis visibility.
- Low profile delivery system makes it potentially easier to reach and treat challenging lesions.

---

**The GORE® VIABAHN®  
Device has a reported  
fracture rate of < .015%  
across all uses.**

(Data on file 2022;  
W. L. Gore & Associates, Inc.; Flagstaff, AZ.)

---

\* CBAS® Heparin Surface. W. L. Gore & Associates Web site. <https://www.goremedical.com/cbas/references>. Accessed May 20, 2019.

# Lasting thromboresistance. Proven technology.\*

The CBAS® Heparin Surface of the GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface\* consists of a proprietary covalent end-point bond that preserves the active site, thus retaining heparin's anticoagulant activity.



#### Proven heparin availability

Performance-ready heparin active site.<sup>21,22</sup>



#### Proven heparin bioactivity

Unmatched, persistent ability to take up antithrombin<sup>23,24</sup>



#### Proven lasting thromboresistance

Improved surface hemocompatibility resulting from heparin availability and bioactivity.<sup>21-25</sup>

## GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface



The bioactive luminal surface of a 5 mm diameter GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface appears free of thrombus after two hours in an in vitro blood loop model.

## Control endoprosthesis



The non-bioactive luminal surface of a control endoprosthesis (5 mm diameter) appears covered with thrombus after 90 minutes in the same blood loop model (Data on file; W. L. Gore & Associates, Inc; Flagstaff, AZ).

\* CBAS® Heparin Surface. W. L. Gore & Associates Web site. <https://www.goremedical.com/eu/cbas/references>. Accessed September 25, 2019.

# Sizing tables

## GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface

### 0.035" Guidewire compatible

Device sizing		Introducer Sheath (Fr)						
Endoprosthesis labeled diameter (mm)*	Recommended vessel diameter† (mm)	2.5 cm device length*	5 cm device length*	7.5 cm device length*	10 cm device length*	15 cm device length*	25 cm device length*	Recommended balloon diameter for device touch-up‡ (mm)
5	4.0–4.7	7	7	7	7	7	7	5
6	4.8–5.5	7	7	7	7	7	7	6
7	5.6–6.5	8	8	8	8	8	8	7
8	6.6–7.5	8	8	8	8	8	8	8
9	7.6–8.5	–	8	8	8	8	–	9
10	8.6–9.5	–	8	–	8	8	–	10
11	9.6–10.5	–	10	–	10	–	–	12
13	10.6–12.0	–	10‡	–	10‡	–	–	14

### 0.018" Guidewire compatible

Device sizing		Introducer Sheath (Fr)						
Endoprosthesis labeled diameter (mm)*	Recommended vessel diameter† (mm)	2.5 cm device length*	5 cm device length*	7.5 cm device length*	10 cm device length*	15 cm device length*	25 cm device length*	Recommended balloon diameter for device touch-up‡ (mm)
5	4.0–4.7	6	6	6	6	6	6	5
6	4.8–5.5	6	6	6	6	6	6	5
7	5.6–6.5	7	7	7	7	7	7	7
8	6.6–7.5	7	7	7	7	7	7 <sup>¶</sup>	8

\* Labeled device diameters and lengths are nominal.

† Recommended endoprosthesis compression within the vessel is approximately 5 – 20%.

‡ The 13 mm diameter device is not compatible with the 10 Fr COOK® FLEXOR® CHECK-FLO® Sheath.

§ For the 11 and 13 mm diameter devices, balloon inflation pressure should not exceed 8 atm.

¶ The 8 mm x 25 cm device is not compatible with the 7 Fr COOK® FLEXOR® CHECK-FLO® Sheath.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

## Advancing Care Through Access



A new standard of flexibility to treat more challenging anatomies

- Deliver with ease: Hydrophilic coating and enhanced flexibility provide exceptional access to challenging anatomies and branch vessels
- Minimize blood loss: Exclusive GORE® DrySeal valve enables introduction of multiple devices with proven hemostasis control

### GORE® DRYSEAL Flex Introducer Sheath — 10 Fr

Catalogue number	Sheath size (Fr)	Minimum sheath ID (mm)	Nominal sheath OD (mm)	Working length (cm)
DSF1033	10	3.3	4.0	33
DSF1045	10	3.3	4.0	45
DSF1065	10	3.3	4.0	65

## References

1. Piazza M, Squizzato F, Dall'Antonia A, et al., Outcomes of self expanding PTFE covered stent versus bare metal stent for chronic iliac artery occlusion in matched cohorts using propensity score modelling. *European Journal of Vascular & Endovascular Surgery* 2017;54(2):177-185.
2. Chang RW, Goodney PP, Baek JH, Nolan BW, Rzucidlo EM, Powell RJ. Long-term results of combined common femoral endarterectomy and iliac stenting/stent grafting for occlusive disease. *Journal of Vascular Surgery* 2008;48(2):362-367.
3. Vesely T, DaVanzo W, Behrend T, Dwyer A, Aruny J. Balloon angioplasty versus Viabahn stent graft for treatment of failing or thrombosed prosthetic hemodialysis grafts. *Journal of Vascular Surgery* 2016;64(5):1400-1410.e1. <http://www.sciencedirect.com/science/article/pii/S0741521416301756>
4. Vesely T, Rodríguez A. Summary of the Gore REVISE Clinical Study. *Endovascular Today* 2014;13(6)Supplement:22-26.
5. Mohr BA, Sheen AL, Roy-Chaudhury P, Schultz SR, Aruny JE; REVISE Investigators. Clinical and economic benefits of stent grafts in dysfunctional and thrombosed hemodialysis access graft circuits in the REVISE Randomized Trial. *Journal of Vascular & Interventional Radiology* 2018;30(2):203-211.e4.
6. W. L. Gore & Associates, Inc. GORE® VIABAHN® Endoprosthesis versus Percutaneous Transluminal Angioplasty (PTA) to Revise Arteriovenous Grafts at the Venous Anastomosis in Hemodialysis Patients. (GORE REVISE Study, AVR 06-01). Flagstaff, AZ: W. L. Gore & Associates, Inc; 2012. [IDE Final Clinical Study Report]. G070069.
7. Ohki T, Kichikawa K, Yokoi H, et al., Outcomes of the Japanese multicenter Viabahn trial of endovascular stent grafting for superficial femoral artery lesions. *Journal of Vascular Surgery* 2017;66(1):130-142.e1.
8. Lammer J, Zeller T, Hausegger KA, et al., Sustained benefit at 2 years for covered stents versus bare-metal stents in long SFA lesions: the VIASTAR Trial. *Cardiovascular & Interventional Radiology* 2015;38(1):25-32.
9. Zeller T, Peeters P, Bosiers M, et al., Heparin-bonded stent-graft for the treatment of TASC II C and D femoropopliteal lesions: the Viabahn-25 cm Trial. *Journal of Endovascular Therapy* 2014;21(6):765-774.
10. Reijnen MMPJ, van Walraven LA, Fritschy WM, et al., 1-year results of a multicenter randomized controlled trial comparing heparin-bonded endoluminal to femoropopliteal bypass. *JACC: Cardiovascular Interventions* 2017;10(22):2320-2331. <http://www.sciencedirect.com/science/article/pii/S1936879817319775>.
11. Saxon RR, Chervu A, Jones PA, et al., Heparin bonded, expanded polytetrafluoroethylene lined stent graft in the treatment of femoropopliteal artery disease: 1 year results of the VIPER (Viabahn Endoprosthesis with Heparin Bioactive Surface in the Treatment of Superficial Femoral Artery Obstructive Disease) Trial. *Journal of Vascular & Interventional Radiology* 2013;24(2):165-173.
12. Yamaoka T. VIABAHN the latest real world clinical data from Japan to the worlds PMS 1Y/IDE 5Y VIABAHN. Presented at JETTALKS on Air; April 18-19, 25-26, 2020; Osaka, Japan.
13. Iida O, Takahara M, Soga Y, et al; VANQUISH Investigators. One-year outcomes of heparin-bonded stent-graft therapy for real-world femoropopliteal lesions and the association of patency with the prothrombotic state based on the prospective, observational, multicenter Viabahn Stent-Graft Placement for Femoropopliteal Disease Requiring Endovascular Therapy (VANQUISH) Study. *Journal of Endovascular Therapy* 2021;28(1):123-131.
14. McQuade K, Gable D, Pearl G, Theune B, Black S. Four-year randomized prospective comparison of percutaneous ePTFE/nitinol self-expanding stent graft versus prosthetic femoral-popliteal bypass in the treatment of superficial femoral artery occlusive disease. *Journal of Vascular Surgery* 2010;52(3):584-591.
15. GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface Instructions for Use (IFU). W. L. Gore & Associates, Inc website. Accessed July 28, 2020. <https://eifu.goremedical.com/7>.
16. Iida O. Update on value and indications of the Viabahn self-expanding stent-graft for fempop occlusive disease: evolution of the device: technical tips and 5-year results from Japan. Presented at the 46th Annual Symposium on Vascular and Endovascular Issues, Techniques, Horizons (VEITHsymposium); November 19 -23, 2019; New York, NY.
17. Troutman DA, Madden NJ, Dougherty MJ, Calligaro KD. Duplex ultrasound diagnosis of failing stent grafts placed for occlusive disease. *Journal of Vascular Surgery* 2014;60(6):1580-1584.
18. Bosiers M, DeLoose K, Callaert J, et al., Superiority of stent-grafts for in-stent restenosis in the superficial femoral artery: twelve-month results from a multicenter randomized trial. *Journal of Endovascular Therapy* 2015;22(1):1-10.
19. Golchehr B, Zeebregts CJ, Reijnen MMPJ, Tielliu IFJ. Long-term outcome of endovascular popliteal artery aneurysm repair. *Journal of Vascular Surgery*. 2018; 67(6): 1797-1804.
20. Beuschel B, Nayfeh T, Kunbaz A, et al., A systematic review and meta-analysis of treatment and natural history of popliteal artery aneurysms. *Journal of Vascular Surgery* 2022;75(1)Supplement:1215-1255.e14.
21. Gore S, Andersson J, Biran R, Underwood C, Riesenfeld J. Heparin surfaces: impact of immobilization chemistry on hemocompatibility and protein adsorption. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 2014;102(8):1817-1824.
22. Biran R, Pond D. Heparin coatings for improving blood compatibility of medical devices. *Advanced Drug Delivery Reviews* 2017;112:12-23.
23. Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® BioActive Surface heparin immobilization. *European Journal of Vascular & Endovascular Surgery* 2003;25(5):432-437.
24. Freeman J, Chen A, Weinberg RJ, Okada T, Chen C, Lin PH. Sustained thromboresistant bioactivity with reduced intimal hyperplasia of heparin-bonded PTFE Propaten Graft in a chronic canine femoral artery bypass model. *Annals of Vascular Surgery* 2018;49:295-303.
25. CBAS® Heparin Surface. W. L. Gore & Associates Web site. <https://www.goremedical.com/eu/cbas/references>. Accessed September 25, 2019.

 Consult Instructions  
for Use  
**eifu.goremedical.com**

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

GORE, *Together, improving life*, DRYSEAL, HEMOBAHN, PROPATEN, VIABAHN and designs are trademarks of W. L. Gore & Associates.

© 2023 W. L. Gore & Associates GmbH 231220052-EN OCTOBER 2023

**W. L. Gore & Associates, Inc.**  
goremedical.com

**Asia Pacific** +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673  
**United States** Flagstaff, AZ 86003 800 437 8181 928 779 2771

