GORE® SYNECOR

Intraperitoneales Biomaterial

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
Linn JG, Mallico EJ, Doerhoff CR, Grantham DW, Washington RG Jr. Evaluation of long-term performance of an intraperitoneal biomaterial in the treatment of ventral hernias. *Surgical Endoscopy* 2023;37:3455–3462. <https://link.springer.com/article/10.1007/s00464-022-09803-9>

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Literatur

1. Linn JG, Mallico EJ, Doerhoff CR, Grantham DW, Washington RG Jr. Evaluation of long-term performance of an intraperitoneal biomaterial in the treatment of ventral hernias. *Surgical Endoscopy* 2023;37:3455-3462.
2. W. L. Gore & Associates, Inc. *Clinical Evaluation Report for GORE® SYNECOR Intraperitoneal Biomaterial*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2022. [Klinischer Bewertungsbericht (Clinical Evaluation Report, CER)]. MD146932. 5. Fassung.
3. Petro CC, Nahabet EH, Criss CN, et al. Central failures of lightweight monofilament polyester mesh causing hernia recurrence: a cautionary note. *Hernia* 2015;19(1):155-159.
4. Cobb WS, Warren JA, Ewing JA, Burnikel A, Merchant M, Carbonell AM. Open retromuscular mesh repair of complex incisional hernia: predictors of wound events and recurrence. *Journal of the American College of Surgeons* 2015;220(4):606-613.
5. Warren JA, McGrath SP, Hale AL, Ewing JA, Carbonell AM 2nd, Cobb WS 4th. Patterns of recurrence and mechanisms of failure after open ventral hernia repair with mesh. *American Surgeon* 2017;83(11):1275-1282.
6. W. L. Gore & Associates, Inc. *GORE SYNECOR Intraperitoneal Biomaterial Design Control (DC) Matrix*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Design Control Matrix – DC Matrix]. MD184278. 2. Fassung.
7. Ventral Hernia Working Group, Breuing K, Butler CE, et al. Incisional ventral hernias: review of the literature and recommendations regarding the grading and technique of repair. *Surgery* 2010;148(3):544-558.
8. Crawford N. *Assessment of Vascularity via Micro CT in Various Patch Devices*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Abschließender Studienbericht]. 2344TL.
9. Berman A. *Evaluation of Plexus with film and ETHICON PHYSIOMESH® in a 30-day rabbit intraperitoneal model*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2015. [Studienprotokoll]. 23355C.
10. W. L. Gore & Associates, Inc. *Plexus Knit PQ Validation Report*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2022. [Validierungsbericht]. MD145325. 5. Fassung.
11. Olson TB. *Competitive Hernia Device Strength Evaluation*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Arbeitsplan]. WP108484.
12. Matthews BD. Absorbable and nonabsorbable barriers on prosthetic biomaterials for adhesion prevention after intraperitoneal placement of mesh. *International Surgery* 2005;90(3)Supplement:530-534.
13. Berman A. *Evaluation of Plexus with film and ETHICON PHYSIOMESH® in a 180-day rabbit intraperitoneal model*. Flagstaff, AZ: W.L. Gore & Associates, Inc; 2015. [Studienprotokoll]. 23375C.
14. Sharkawy AA, Klitzman B, Truskey GA, Reichert WM. Engineering the tissue which encapsulates subcutaneous implants. II. Plasma-tissue exchange properties. *Journal of Biomedical Materials Research* 1998;40(4):586-597.
15. Holt DJ, Grainger DW. Host response to biomaterials. In: Hollinger JO, ed. *An Introduction to Biomaterials*. 2nd ed. Boca Raton, FL: CRC Press; 2012;6:91-118.
16. Berman A. *A Rabbit Model for the Biomechanical and Histological Assessment of Suture Line Wound Healing*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2008. [Studienprotokoll]. 19785C.
17. Morales-Conde S, Flores M, Fernández V, Morales-Méndez S. Bioabsorbable vs polypropylene plug for the “Mesh and Plug” inguinal hernia repair. Poster präsentiert beim 9. Jahrestreffen der American Hernia Society. 9.-12. Februar 2005, San Diego, CA.
18. LeBlanc KA, Heniford BT, Voeller GR. Innovations in ventral hernia repair. Materials and techniques to reduce MRSA and other infections. *Contemporary Surgery* 2006;62(4)Supplement:1-8.
19. Clinger L. *PTFE Knit Microbial Placement*. Flagstaff, AZ; W. L. Gore & Associates, Inc; 2018. [Arbeitsplan]. WP110158.
20. Lake SP, Ray S, Zihni AM, Thompson DM Jr, Gluckstein J, Deeken CR. Pore size and pore shape—but not mesh density—alter the mechanical strength of tissue ingrowth and host tissue response to synthetic mesh materials in a porcine model of ventral hernia repair. *Journal of the Mechanical Behavior of Biomedical Materials* 2015; 42:186-197.
21. Klinge U, Klosterhalfen B, Birkenhauer V, Junge K, Conze J, Schumpelick V. Impact of polymer pore size on the interface scar formation in a rat model. *Journal of Surgical Research* 2002;103(2):208-214.
22. Zhu LM, Schuster P, Klinge U. Mesh implants: an overview of crucial mesh parameters. *World Journal of Gastrointestinal Surgery* 2015;7(10):226-236.
23. Klinge U, Klosterhalfen B, Conze J, et al. Modified mesh for hernia repair that is adapted to the physiology of the abdominal wall. *European Journal of Surgery* 1998;164(12):951-960.
24. Socea B, Socea LI, Bratu OG, et al. Recurrence rates and mesh shrinkage after polypropylene vs. polyester mesh hernia repair in complicated hernias. *Materiale Plastice* 2018;55(1):79-81.
25. Klosterhalfen B, Junge K, Klinge U. The lightweight and large porous mesh concept for hernia repair. *Expert Reviews in Medical Devices* 2005;2(1):103-117.
26. Klitchko M. *SYNECOR IP - Comparison of Film Strength to Competitive Intraperitoneal Hernia Meshes*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2019. [Arbeitsplan]. WP110536.
27. Jenkins ED, Yip M, Melman L, Frisella MM, Matthews BD. Informed consent: cultural and religious issues associated with the use of allogeneic and xenogeneic mesh products. *Journal of the American College of Surgeons* 2010;210(4):402-410.

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© 2023, 2024 W. L. Gore & Associates GmbH 241354228-DE JANUAR 2024
Inc. TELA BIO und OVITEX sind Markenzeichen von TELA Bio, Inc.