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
Daenens K, Schepers S, Fourneau I, Houthoofd S, Nevelsteen A. Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results. *Journal of Vascular Surgery* 2009;49(5):1210-1216.

Hugl B, Nevelsteen A, Daenens K, *et al*; PEPE II Study Group. PEPE II - a multicenter study with an end-point heparin-bonded expanded polytetrafluoroethylene vascular graft for above and below knee bypass surgery: determinants of patency. *Journal of Cardiovascular Surgery* 2009;50(2):195-203.

Kirkwood ML, Wang GJ, Jackson BM, Golden MA, Fairman RM, Woo EY. Lower limb revascularization for PAD using a heparin-coated PTFE conduit. *Vascular & Endovascular Surgery* 2011;45(4):329-334.

Lösel-Sadée H, Alefelder C. Heparin-bonded expanded polytetrafluoroethylene graft for infragenicular bypass: 5-year results. *Journal of Cardiovascular Surgery* 2009;50(3):339-343.

Peeters P, Verbist J, Deloose K, Bosiers M. Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass? *Italian Journal of Vascular & Endovascular Surgery* 2008;15(3):143-148.

Pulli R, Dorigo W, Castelli P, *et al*; Propaten Italian Registry Group. Midterm results from a multicenter registry on the treatment of infrainguinal critical limb ischemia using a heparin-bonded ePTFE graft. *Journal of Vascular Surgery* 2010;51(5):1167-1177.

Walluscheck KP, Bierkandt S, Brandt M, Cremer J. Infrainguinal ePTFE vascular graft with bioactive surface heparin bonding-first clinical results. *Journal of Cardiovascular Surgery* 2005;46(4):425-430.

Battaglia G, Tringale R, Monaca V. Retrospective comparison of a heparin bonded ePTFE graft and saphenous vein for infragenicular bypass: implications for standard treatment protocol. *Journal of Cardiovascular Surgery* 2006;47(1):41-47.

Bosiers M, Deloose K, Verbist J, *et al*. Heparin-bonded expanded polytetrafluoroethylene vascular graft for femoropopliteal and femorocrural bypass grafting: 1-year results. *Journal of Vascular Surgery* 2006;43(2):313-319.

Dorigo W, Di Carlo F, Troisi N, *et al*. Lower limb revascularization with a new bioactive prosthetic graft: early and late results. *Annals of Vascular Surgery* 2008;22(1):79-87.

Dorigo W, Pulli R, Alessi Innocenti A, *et al*. Lower limb below-knee revascularization with a new bioactive prosthetic graft: a case-control study. *Italian Journal of Vascular & Endovascular Surgery* 2005;12(3):75-81.

Dorigo W, Pulli R, Castelli P, *et al*; Propaten Italian Registry Group. A multicentric comparison between autologous saphenous vein and heparin-bonded expanded polytetrafluoroethylene (ePTFE) graft in the treatment of critical limb ischemia in diabetics. *Journal of Vascular Surgery* 2011;54(5):1332-1338.

Dorrucci V, Griselli F, Petralia G, Spinamano L, Adornetto R. Heparin-bonded expanded polytetrafluoroethylene grafts for infragenicular bypass in patients with critical limb ischemia: 2 year results. *Journal of Cardiovascular Surgery* 2008;49(2):145-149.

Lindholt JS, Gottschalksen B, Johannesen N, *et al*. The Scandinavian Propaten® Trial – 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses – a randomised clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673.

Peeters P, Verbist J, Deloose K, Bosiers M. Results with heparin bonded PTFE grafts for femorodistal bypasses. *Journal of Cardiovascular Surgery* 2006;47(4):407-413.

## Referências

1. Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® BioActive Surface heparin immobilization. *European Journal of Vascular & Endovascular Surgery* 2003;25(5):432-437.
2. Fisher JL, Thomson RC, Moore JW, Begovac PC. Functional parameters of thromboresistant heparinized e-PTFE vascular grafts. *Cardiovascular Pathology* 2002;11(1):42.
3. Heyligers JMM, Lisman T, Weeterings C, et al. Heparin immobilization reduces thrombogenicity on small-caliber ePTFE grafts. *Journal of Vascular Surgery* 2006;43:587-591.
4. Lin PH, Chen C, Bush RL, Yao Q, Lumsden AB, Hanson SR. Small-caliber heparin-coated ePTFE grafts reduce platelet deposition and neointimal hyperplasia in a baboon model. *Journal of Vascular Surgery* 2004; 39(6):1322-1328.
5. Lin PH, Bush RL, Yao Q, Lumsden AB, Chen C. Evaluation of platelet deposition and neointimal hyperplasia of heparin-coated small-caliber ePTFE grafts in a canine femoral artery bypass model. *Journal of Surgical Research* 2004;118:45-52.
6. Kocsis JF, Llanos G, Holmer E. Heparin-coated stents. *Journal of Long-Term Effects of Medical Implants* 2000;10(1&2):19-45.
7. Chen C, Lumsden AB, Hanson SR. Local infusion of heparin reduces anastomotic neointimal hyperplasia in aortoiliac expanded polytetrafluoroethylene bypass grafts in baboons. *Journal of Vascular Surgery* 2000;31(2):354-363.
8. Belch JJ, Dormandy J, Biasi GM, et al; CASPAR Writing Committee. Results of the randomized, placebo-controlled clopidogrel and acetylsalicylic acid in bypass surgery for peripheral arterial disease (CASPAR) trial. *Journal of Vascular Surgery* 2010;52(4):825-833.
9. Heyligers JMM, Lisman T, Verhagen HJM, Weeterings C, de Groot PG, Moll FL. A heparin-bonded vascular graft generates no systemic effect on markers of hemostasis activation or detectable heparin-induced thrombocytopenia-associated antibodies in humans. *Journal of Vascular Surgery* 2008;47(2):324-329.
10. Lindholt JS, Gottschalksen B, Johannesen N, et al. The Scandinavian Propaten® Trial – 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses – a randomised clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673.
11. Warkentin TE. An overview of the heparin-induced thrombocytopenia syndrome. *Seminars in Thrombosis & Hemostasis* 2004;30(3):273-283.
12. Franchini M. Heparin-induced thrombocytopenia: an update. *Thrombosis Journal* 2005;3(1):14-18.
13. Jang IK, Hursting MJ. When heparins promote thrombosis: review of heparin-induced thrombocytopenia. *Circulation* 2005;111(20):2671-2683.
14. Chong BH. Heparin-induced thrombocytopenia. *Journal of Thrombosis & Haemostasis* 2003;1(7):1471-1478.
15. Lubenow N, Kempf R, Eichner A, Eichler P, Carlsson LE, Greinacher A. Heparin-induced thrombocytopenia: temporal pattern of thrombocytopenia in relation to initial use or reexposure to heparin. *Chest* 2002;122(1):37-42.
16. Warkentin TE, Kelton JG. Temporal aspects of heparin-induced thrombocytopenia. *New England Journal of Medicine* 2001;344(17):1286-1292.
17. Jackson MR, Neilson WJ, Lary M, Baay P, Web K, Clagett GP. Delayed-onset heparin-induced thrombocytopenia and thrombosis after intraoperative heparin anticoagulation-four case reports. *Vascular & Endovascular Surgery* 2006;40(1):67-70.
18. Rice L, Attisha WK, Drexler A, Francis JL. Delayed-onset heparin-induced thrombocytopenia. *Annals of Internal Medicine* 2002;136(3):210-215.
19. Smythe MA, Stephens JL, Mattson JC. Delayed-onset heparin-induced thrombocytopenia. *Annals of Emergency Medicine* 2005;45(4):417-419.
20. Warkentin TE, Kelton JG. Delayed-onset heparin-induced thrombocytopenia and thrombosis. *Annals of Internal Medicine* 2001;135(7):502-506.
21. Kasirajan K. Outcomes following heparin-induced thrombocytopenia in patients with heparin bonded vascular grafts. Abstract presented at the 2011 Vascular Annual Meeting; June 16-18, 2011; Chicago, IL. *Journal of Vascular Surgery* 2011;53(6)Supplement 1:17S-18S. Abstract SS4.
22. Warkentin TE, Greinacher A, Koster A, Lincoff AM; American College of Chest Physicians. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines (8th Edition). *Chest* 2008;133(6)Supplement:340S-380S.

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GORE®, GORE-TEX®, PERFORMANCE THROUGH COLLABORATION, PROPATEN®, and designs are trademarks of W. L. Gore & Associates.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

CORDIS®, BX VELOCITY®, and HEPACOAT® are trademarks of Cordis, a Johnson & Johnson Company. Medline® is a registered trademark of the U.S. National Library of Medicine. All other trademarks / tradenames are property of their respective owners.

© 2014, 2016 W. L. Gore & Associates, Inc. AK0409-PT4 JULY 2016