

## THE VANQUISH STUDY<sup>1</sup>


A prospective, multicenter, real-world study of the VIABAHN® Device in femoropopliteal occlusive disease

Large real-world study of the VIABAHN® Device

**343**

Lesions in the full coverage cohort

Proven patency



Consistent with other endovascular therapies, sole risk factor for loss of patency was small vessel diameter (as measured by intravascular ultrasound (IVUS)).

Challenging disease

**25 cm**

Mean lesion length

**71.4%**

Chronic total occlusions (CTO)

**23.3%**

Chronic limb ischemia

**90.1%**

TASC II C (39.4%) and D (50.7%)

**80% average primary patency at one year in over 1,000 lesions<sup>1-7</sup>**

Proven patency in complex superficial femoral artery lesions across seven multicenter, prospective, randomized or single-arm studies

Gore VIPER Clinical Study <sup>2</sup> N = 119	VIASTAR Trial <sup>3</sup> N = 66	25 cm Trial <sup>4</sup> N = 71	SuperB Study <sup>5</sup> N = 63	Gore Japan IDE Clinical Study <sup>6</sup> N = 103	Gore Japan Post Market Clinical Study <sup>7</sup> N = 324	<b>VANQUISH Study<sup>1</sup></b> N = 343
<b>19 cm</b> Mean lesion length	<b>19 cm</b> Mean lesion length	<b>27 cm</b> Mean lesion length	<b>23 cm</b> Mean lesion length	<b>22 cm</b> Mean lesion length	<b>24 cm</b> Mean lesion length	<b>25 cm</b> Mean lesion length
<b>56%</b> CTO	<b>79%</b> CTO	<b>93%</b> CTO	<b>75%</b> CTO <sup>†</sup>	<b>66%</b> CTO	<b>70%</b> CTO <sup>†</sup>	<b>71%</b> CTO

\* As used by Gore, Heparin Bioactive Surface refers to Gore's proprietary CBAS Heparin Surface.

## References and citations

† CTO percentage defined as percentage of TASC II D.

‡ Data on file 2020; W. L. Gore & Associates, Inc; Flagstaff AZ.

1. Iida O, Takahara M, Soga Y, et al; VANQUISH Investigators. One-year outcomes of heparin-bonded stent-graft therapy for real-world femoropopliteal lesions and the association of patency with the prothrombotic state based on the prospective, observational, multicenter Viabahn Stent-Graft Placement for Femoropopliteal Diseases Requiring Endovascular Therapy (VANQUISH) Study. *Journal of Endovascular Therapy* 2021;28(1):123-131.
2. Saxon RR, Chervu A, Jones PA, et al. Heparin bonded, expanded polytetrafluoroethylene lined stent graft in the treatment of femoropopliteal artery disease: 1 year results of the VIPER (Viabahn Endoprosthesis with Heparin Bioactive Surface in the Treatment of Superficial Femoral Artery Obstructive Disease) Trial. *Journal of Vascular & Interventional Radiology* 2013;24(2):165-173.
3. Lammer J, Zeller T, Hausegger KA, et al. Sustained benefit at 2 years for covered stents versus bare-metal stents in long SFA lesions: the VIASTAR Trial. *Cardiovascular & Interventional Radiology* 2015;38(1):25-32.
4. Zeller T, Peeters P, Bosiers M, et al. Heparin-bonded stent-graft for the treatment of TASC II C and D femoropopliteal lesions: the Viabahn-25 cm Trial. *Journal of Endovascular Therapy* 2014;21(6):765-774.
5. Reijnen MMPJ, van Walraven LA, Fritschy WM, et al. 1-year results of a multicenter randomized controlled trial comparing heparin-bonded endoluminal to femoropopliteal bypass. *JACC: Cardiovascular Interventions* 2017;10(22):2320-2331.
6. Ohki T, Kichikawa K, Yokoi H, et al. Outcomes of the Japanese multicenter Viabahn trial of endovascular stent grafting for superficial femoral artery lesions. *Journal of Vascular Surgery* 2017;66(1):130-142.e1.
7. Yamaoka T. VIABAHN the latest real world clinical data from Japan to the worlds PMS 1Y/IDE 5Y VIABAHN. Presented at JETTALKS on Air; April 18-19, 25-26, 2020; Osaka, Japan.

 Consult Instructions  
for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

**INDICATIONS FOR USE IN THE U.S.:** The GORE® VIABAHN® Endoprosthesis is indicated for improving blood flow in patients with symptomatic peripheral arterial disease in superficial femoral artery de novo and restenotic lesions up to 270 mm in length with reference vessel diameters ranging from 4.0 – 7.5 mm. The GORE® VIABAHN® Endoprosthesis is indicated for improving blood flow in patients with symptomatic peripheral arterial disease in superficial femoral artery in-stent restenotic lesions up to 270 mm in length with reference vessel diameters ranging from 4.0 – 6.5 mm. The GORE® VIABAHN® Endoprosthesis is indicated for improving blood flow in patients with symptomatic peripheral arterial disease in iliac artery lesions up to 80 mm in length with reference vessel diameters ranging from 4.0 – 12 mm. The GORE® VIABAHN® Endoprosthesis is also indicated for the treatment of stenosis or thrombotic occlusion at the venous anastomosis of synthetic arteriovenous (AV) access grafts. **CONTRAINDICATIONS:** The GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface is contraindicated for non-compliant lesions where full expansion of an angioplasty balloon catheter was not achieved during pre-dilatation, or where lesions cannot be dilated sufficiently to allow passage of the delivery system. Do not use the GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface in patients with known hypersensitivity to heparin, including those patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II. Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available.  $\text{Rx}_{\text{Only}}$

Products listed may not be available in all markets.

CBAS is a trademark of Carmeda, AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.  
GORE, *Together, improving life*, VIABAHN and designs are trademarks of W. L. Gore & Associates.  
© 2021 W. L. Gore & Associates, Inc. 2019385-EN JULY 2021

W. L. Gore & Associates, Inc.  
[goremedical.com](http://goremedical.com)

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673  
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

