




GORE® ACUSEAL
Vascular Graft

UNCOMPROMISED HANDLING WITH TRI-LAYER SEALING PROPERTIES

Cannulation capable
within 24 hours



Together, improving life


Low bleed barrier

- Elastomeric middle layer
- Low-bleed through puncture sites, hinders cannulation needle bleeding
- Hinders suture line bleeding
- May reduce risk of seroma formation*

Low bleed versus bleed




Low bleed: GORE® ACUSEAL Vascular Graft



Bleed: Standard ePTFE Graft

Post cannulation of the luminal surface with a 16 gauge needle. Hold pressure for 10–15 minutes to achieve hemostasis post needle removal.

Tri-layer construction of a GORE® ACUSEAL Vascular Graft



Abluminal layer: ePTFE graft

Elastomeric layer

Luminal layer: ePTFE with CBAS Heparin Surface


* Data on file.

Uncompromised handling

- Flexible at curves without kinking
- Free from stiffness or rigidity
- Precise suturing and anastomotic tailoring




GORE® ACUSEAL Vascular Graft with cannulation needle through graft wall.




GORE® ACUSEAL Vascular Graft: flexibility without kinking.

A thromboresistant luminal graft surface

Evaluation of GORE® ACUSEAL Vascular Graft in a benchtop canine blood flow loop model.



GORE® ACUSEAL Vascular Graft with CBAS Heparin Surface



GORE® ACUSEAL Vascular Graft without CBAS Heparin Surface

Cannulation capable within 24 hours

- Tri-layer design is optimized for early cannulation
- Expands treatment options for earlier removal or avoidance of a central venous catheter

GORE® ACUSEAL Vascular Graft Clinical Study results* (N = 138)

Cumulative patency	GORE® ACUSEAL Vascular Graft	Historical control
6-month follow-up	84%	75%
12-month follow-up	78%	66%

54 patients (40%) were cannulated within 72 hours of implantation


Time from implantation to first cannulation	Number of GORE® ACUSEAL Vascular Grafts cannulated†
Within 24 hours	n = 30 (22.2%)
Within 48 hours	n = 48 (35.6%)
Within 72 hours	n = 54 (40.0%)
Within 7 days	n = 70 (51.9%)

* Data on file; Flagstaff, Arizona, 2013.

† N = 138, three grafts were not cannulated.


Proprietary covalent end-point bonding

Covalent end-point bonding allows the heparin to extend into the bloodstream, keeping the active site bioavailable, unlike a non-permanent bond that can be washed away in the bloodstream.



- The anticoagulant function of heparin is dependent on the bioavailability of an active site within the molecule.
- Some methods of covalent heparin bonding damage and/or obstruct the active site, and hence destroy heparin's anticoagulant activity.
- The CBAS Heparin Surface of the GORE® ACUSEAL Vascular Graft consists of a proprietary covalent end-point bond that preserves the active site, thus retaining heparin's anticoagulant activity.

Mechanism of action



- Bioactive site of the heparin molecule enables antithrombin to bind thrombin.
- When antithrombin binds to thrombin, a neutral AT-T complex is formed.
- Neutral AT-T complex detaches from the heparin molecule. Active site becomes available to again bind antithrombin.

Products listed may not be available in all markets.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.
GORE, *Together, improving life*, ACUSEAL and designs are trademarks of W. L. Gore & Associates.
© 2021 W. L. Gore & Associates, Inc. 21259650-EN OCTOBER 2021

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

