

GORE® PROPATEN® Vascular Graft

Proven patency. Measurable value.

LITERATURE SUMMARY




Table of contents

| | | | |
|----------------------------------|------|---|-----------|
| Below-knee bypass summary | | | 3 |
| Kaisar et al | 2018 | Comparison of propaten heparin-bonded vascular graft with distal anastomotic patch versus autogenous saphenous vein graft in tibial artery bypass | 4 |
| Uhl et al | 2015 | Comparison of venous and HePTFE tibial and peroneal bypasses in critical limb ischemia patients unsuitable for endovascular revascularization | 5 |
| Neville et al | 2014 | Heparin-bonded ePTFE (Propaten): is it as good as autologous vein for tibial bypass? | 6 |
| Monaca et al | 2013 | Subpopliteal revascularization. Criteria analysis for use of E-PTFE (Propaten®) as first choice conduit | 7 |
| Dorigo et al | 2012 | Results from an Italian multicentric registry comparing heparin-bonded ePTFE graft and autologous saphenous vein in below-knee femoro-popliteal bypasses | 8 |
| Pulli et al | 2010 | Midterm results from a multicenter registry on the treatment of infrainguinal critical limb ischemia using a heparin-bonded ePTFE graft | 9 |
| Daenens et al | 2009 | Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results | 10 |
| Hugl et al | 2009 | PEPE II — A multicenter study with an end-point heparin-bonded expanded polytetrafluoroethylene vascular graft for above and below knee bypass surgery: determinants of patency | 11 |
| Lösel-Sadée & Alefelder | 2009 | Heparin-bonded expanded polytetrafluoroethylene graft for infragenicular bypass: five-year results | 12 |
| Peeters et al | 2008 | Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass? | 13 |
| Above-knee bypass summary | | | 14 |
| Piffaretti et al | 2018 | Results from a multicenter registry of heparin-bonded expanded polytetrafluoroethylene graft for above-the-knee femoropopliteal bypass? | 15 |
| Samson et al | 2016 | Heparin-bonded expanded polytetrafluoroethylene femoropopliteal bypass grafts outperform expanded polytetrafluoroethylene grafts without heparin in a long-term comparison | 16 |
| Lindholt et al | 2011 | The Scandinavian Propaten® Trial — 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses — a randomized clinical controlled multi-centre trial | 17 |
| Daenens et al | 2009 | Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results | 18 |
| Peeters et al | 2008 | Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass? | 19 |
| Bosiers et al | 2006 | Heparin-bonded expanded polytetrafluoroethylene vascular graft for femoropopliteal and femorocrural bypass grafting: 1-year results | 20 |
| Bibliography | | | 21 |


Glossary of terms

| | | | |
|--------------|---|-----------------|---|
| AK | Above-knee bypass | ESRD | End-stage renal disease |
| BK | Below-knee bypass | HePTFE | Heparin-bonded expanded polytetrafluoroethylene |
| CBAS | A trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc., referring to the proven heparin bonding technology on GORE® PROPATEN® Vascular Graft | SFA | Superficial femoral artery |
| CLI | Critical limb ischemia | TP trunk | Tibioperoneal trunk |
| ePTFE | Expanded polytetrafluoroethylene | | |

Legend




Lower limb anatomy



Overall weighted average* primary patency in below-knee bypasses

| | 1 year | 2 years | 3 years | 4 years |
|---|-------------------------------------|---------------------------------------|-------------------------------------|-------------------------------------|
| Vein | 77%[†] N = 2,936 | 58%^{†,‡} N = 1,272 | 63%[†] N = 675 | 62%^{†,‡} N = 615 |
| GORE® PROPATEN® Vascular Graft | 75%[§] N = 971 | 65%[§] N = 782 | 56%[§] N = 703 | 46%[§] N = 693 |
| ePTFE Vascular Graft | 60%[†] N = 2,549 | 46%[†] N = 2,249 | 40%[†] N = 1,941 | 29%[†] N = 599 |



* Weighted Average =
$$\frac{(N_1 \times \text{Primary Patency}_1) + (N_2 \times \text{PP}_2) + \dots + (N_n \times \text{PP}_n)}{N_1 + N_2 + \dots + N_n}$$

† Data based on an analysis of current literature: several MEDLINE® Database searches were performed to identify publications pertaining to ePTFE synthetic vascular graft and vein infragenicular bypasses. Search criteria included (1) articles published from January 2000 to January 2012, (2) key words used were below knee, polytetrafluoroethylene, prosthetic, bypass, patency, (3) articles in English language, (4) N equal or greater than 30 bypasses, (5) clinical publications, (6) reviews, case reports or meta-analysis articles were excluded, (7) articles containing the key word AV access (including synonyms) were excluded. Articles that did not meet the above criteria were deemed ineligible for this analysis. (data on file 2019; W. L. Gore & Associates, Inc; Flagstaff, AZ.)

‡ In studies where 1-year and 3-year patency data were reported, but 2-year patency data were not reported, the 2-year patency rate used in this analysis was interpolated as the average of the 1-year and 3-year patency rates.

§ Below-knee (BK) inclusion criteria for GORE® PROPATEN® Vascular Graft literature used in this analysis were (1) articles in English language, (2) clinical journal articles or book chapters, (3) non-overlapping patient populations, (4) BK bypass primary patency reported for at least 12 months of follow-up and (5) N = 50 or more BK bypasses. Additional exclusion criteria were (1) reviews, case reports or meta-analysis articles and (2) articles containing the key word AV access (including synonyms). (data on file 2019; W. L. Gore & Associates, Inc; Flagstaff, AZ.)

Comparison of propaten heparin-bonded vascular graft with distal anastomotic patch versus autogenous saphenous vein graft in tibial artery bypass¹

Kaisar et al. 2018

Michael E. DeBakey Department of Surgery, Baylor College of Medicine, Houston, TX, USA



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | | | Study details | | | | | | | | | | | | | | | | | | |
|--|----------------------|---------------------|---------------|---|---|--------------|----|----|-----------------|----|----|------------------|----|----|----------|----|----|---|---|----|--|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>1</td> <td>2</td> </tr> <tr> <td>3</td> <td>20</td> <td>32</td> </tr> <tr> <td>4</td> <td>15</td> <td>24</td> </tr> <tr> <td>5</td> <td>18</td> <td>29</td> </tr> <tr> <td>6</td> <td>8</td> <td>13</td> </tr> </tbody> </table> | | | | N | % | 2 | 1 | 2 | 3 | 20 | 32 | 4 | 15 | 24 | 5 | 18 | 29 | 6 | 8 | 13 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center analysis of prospectively collected data All GORE® PROPATEN® Vascular Graft bypasses utilized a distal patch using either autologous vein (58%) or bovine pericardium (42%) All grafts were 6 mm diameter Wound care techniques included aggressive wound debridement, negative-pressure wound closure and intravenous antibiotics Compared to autologous vein, patients receiving a GORE® PROPATEN® Vascular Graft experienced a shorter operative time and length of hospital stay, although the difference did not reach a statistical significance No statistically significant difference in primary patency, secondary patency, or limb salvage between patients receiving autologous vein and GORE® PROPATEN® Vascular Graft at 4 years |
| | N | % | | | | | | | | | | | | | | | | | | | |
| 2 | 1 | 2 | | | | | | | | | | | | | | | | | | | |
| 3 | 20 | 32 | | | | | | | | | | | | | | | | | | | |
| 4 | 15 | 24 | | | | | | | | | | | | | | | | | | | |
| 5 | 18 | 29 | | | | | | | | | | | | | | | | | | | |
| 6 | 8 | 13 | | | | | | | | | | | | | | | | | | | |
| <p>Bypass indication</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>Claudication</td> <td>10</td> <td>16</td> </tr> <tr> <td>Rest pain</td> <td>23</td> <td>37</td> </tr> <tr> <td>Tissue loss</td> <td>29</td> <td>47</td> </tr> </tbody> </table> | | | | N | % | Claudication | 10 | 16 | Rest pain | 23 | 37 | Tissue loss | 29 | 47 | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | |
| Claudication | 10 | 16 | | | | | | | | | | | | | | | | | | | |
| Rest pain | 23 | 37 | | | | | | | | | | | | | | | | | | | |
| Tissue loss | 29 | 47 | | | | | | | | | | | | | | | | | | | |
| Hypertension | Diabetes | Renal failure | | | | | | | | | | | | | | | | | | | |
| 74% N = 46 | 52% N = 32 | 15% N = 9 | | | | | | | | | | | | | | | | | | | |
| <p>Distal anastomotic vessel</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>TP trunk</td> <td>9</td> <td>15</td> </tr> <tr> <td>Anterior tibial</td> <td>14</td> <td>23</td> </tr> <tr> <td>Posterior tibial</td> <td>22</td> <td>36</td> </tr> <tr> <td>Peroneal</td> <td>17</td> <td>27</td> </tr> </tbody> </table> | | | | N | % | TP trunk | 9 | 15 | Anterior tibial | 14 | 23 | Posterior tibial | 22 | 36 | Peroneal | 17 | 27 | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | |
| TP trunk | 9 | 15 | | | | | | | | | | | | | | | | | | | |
| Anterior tibial | 14 | 23 | | | | | | | | | | | | | | | | | | | |
| Posterior tibial | 22 | 36 | | | | | | | | | | | | | | | | | | | |
| Peroneal | 17 | 27 | | | | | | | | | | | | | | | | | | | |

* GORE® PROPATEN® Vascular Graft group.

"Propaten grafts [GORE® PROPATEN® Vascular Graft] with distal anastomotic patch have similar clinical outcomes compared to the saphenous vein graft in tibial artery bypass. Our data support the use of Propaten graft [GORE® PROPATEN® Vascular Graft] with distal anastomotic patch as a viable conduit of choice in patients undergoing tibial artery bypass." — J. Kaisar


Comparison of venous and HePTFE tibial and peroneal bypasses in critical limb ischemia patients unsuitable for endovascular revascularization²

Uhl et al. 2015

Department of Vascular Surgery, Krankenhaus Barmherzige Brüder, Regensburg, Germany



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics | Study details | | | | | | | | | | | | | | |
|--|---------------|----|---|----------|----|----|-----------------|----|----|------------------|----|----|--|----|----|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>4</td> <td>23</td> <td>26</td> </tr> <tr> <td>5</td> <td>66</td> <td>74</td> </tr> <tr> <td>6</td> <td>7</td> <td>8</td> </tr> </tbody> </table> | | N | % | 4 | 23 | 26 | 5 | 66 | 74 | 6 | 7 | 8 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center analysis of prospectively collected data GORE® PROPATEN® Vascular Graft was utilized in patients lacking an adequate autologous vein (diameter less than 3 mm or more than two vein segments required to achieve needed length) 30-day mortality and major amputation rates were both 7% Survival rates at 1 year and 3 years were 77% and 47%, respectively | | |
| | N | % | | | | | | | | | | | | | |
| 4 | 23 | 26 | | | | | | | | | | | | | |
| 5 | 66 | 74 | | | | | | | | | | | | | |
| 6 | 7 | 8 | | | | | | | | | | | | | |
| <p>CLI Hypertension Diabetes</p> <p>100% 90% 46%</p> <p>N = 89 N = 80 N = 41</p> | | | | | | | | | | | | | | | |
| <p>Renal insufficiency Prior operation</p> <p>46% 48%</p> <p>N = 41 N = 43</p> | | | | | | | | | | | | | | | |
| <p>Distal anastomotic vessel</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>TP trunk</td> <td>4</td> <td>5</td> </tr> <tr> <td>Anterior tibial</td> <td>52</td> <td>67</td> </tr> <tr> <td>Posterior tibial</td> <td>12</td> <td>15</td> </tr> <tr> <td>Peroneal</td> <td>10</td> <td>13</td> </tr> </tbody> </table> | | N | % | TP trunk | 4 | 5 | Anterior tibial | 52 | 67 | Posterior tibial | 12 | 15 | Peroneal | 10 | 13 |
| | N | % | | | | | | | | | | | | | |
| TP trunk | 4 | 5 | | | | | | | | | | | | | |
| Anterior tibial | 52 | 67 | | | | | | | | | | | | | |
| Posterior tibial | 12 | 15 | | | | | | | | | | | | | |
| Peroneal | 10 | 13 | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | |

“The results of our study show that autologous vein grafts are still first choice for tibial and peroneal bypasses in patients with critical limb ischemia. If no adequate vein is available, heparin-bonded expanded polytetrafluoroethylene bypasses are an acceptable alternative to an otherwise impending major amputation.” — C. Uhl

Heparin-bonded ePTFE (Propaten): is it as good as autologous vein for tibial bypass?³

Neville et al. 2014

Division of Vascular Surgery, George Washington University, Washington, DC, USA



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | | | Study details | | | | | | | | | | | | | | | |
|---|--------------------------------------|----|---------------|---|---|-----------------|----|----|------------------|----|----|----------|----|----|----------------|----|----|---|
| <table border="1"> <caption>Rutherford classification</caption> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>5</td> <td>8</td> </tr> <tr> <td>4</td> <td>17</td> <td>27</td> </tr> <tr> <td>5</td> <td>30</td> <td>48</td> </tr> <tr> <td>6</td> <td>10</td> <td>16</td> </tr> </tbody> </table> | | | | N | % | 3 | 5 | 8 | 4 | 17 | 27 | 5 | 30 | 48 | 6 | 10 | 16 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center study Comparing GORE® PROPATEN® Vascular Graft to great saphenous vein Distal vein patch was used at all distal anastomoses There was no statistically significant difference in primary patency or limb salvage between GORE® PROPATEN® Vascular Graft and single-segment great saphenous vein |
| | N | % | | | | | | | | | | | | | | | | |
| 3 | 5 | 8 | | | | | | | | | | | | | | | | |
| 4 | 17 | 27 | | | | | | | | | | | | | | | | |
| 5 | 30 | 48 | | | | | | | | | | | | | | | | |
| 6 | 10 | 16 | | | | | | | | | | | | | | | | |
| Hypertension 73% N = 45 | Diabetes 47% N = 29 | | | | | | | | | | | | | | | | | |
| ESRD 13% N = 8 | Prior bypass 44% N = 27 | | | | | | | | | | | | | | | | | |
| <table border="1"> <caption>Distal anastomosis</caption> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>Anterior tibial</td> <td>15</td> <td>24</td> </tr> <tr> <td>Posterior tibial</td> <td>22</td> <td>35</td> </tr> <tr> <td>Peroneal</td> <td>21</td> <td>34</td> </tr> <tr> <td>Dorsalis pedis</td> <td>4</td> <td>6</td> </tr> </tbody> </table> | | | | N | % | Anterior tibial | 15 | 24 | Posterior tibial | 22 | 35 | Peroneal | 21 | 34 | Dorsalis pedis | 4 | 6 | |
| | N | % | | | | | | | | | | | | | | | | |
| Anterior tibial | 15 | 24 | | | | | | | | | | | | | | | | |
| Posterior tibial | 22 | 35 | | | | | | | | | | | | | | | | |
| Peroneal | 21 | 34 | | | | | | | | | | | | | | | | |
| Dorsalis pedis | 4 | 6 | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | |

* GORE® PROPATEN® Vascular Graft group.

“This experience with heparin-bonded ePTFE grafts for solely tibial artery bypass yielded patency and limb salvage rates that are comparable to intact great saphenous vein.” — R. Neville

“We believe that a quality saphenous vein remains the ideal conduit for tibial bypass, although HePTFE should be considered when intact ipsilateral or contralateral vein is not available. In our practice, HePTFE has emerged as the choice over arm vein, especially in the ESRD patient who needs upper extremity vein for dialysis access. We would also choose HePTFE over composite short saphenous vein given the increased dissection required and length of conduit.” — R. Neville

Subpopliteal revascularization. Criteria analysis for the use of E-PTFE (Propaten®) as first choice conduit⁴

Monaca et al. 2013

Vittorio Emanuele Policlinic University Hospital, Presidio Ospedaliero "Ferrarotto", Catania, Italy



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|---------------------------|---------------|--|--|---|---|---|----|----|---|----|----|---|----|----|--------------|----------|---------------|------------|------------|------------|---------|--------|--------|----------------------|--|--|--|---|---|-----------------------|-----|-----|--------------------|--|--|--|---|---|--------------|-----|----|----------|----|----|---|
| <table border="1"> <thead> <tr> <th colspan="3">Rutherford classification</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>87</td> <td>41</td> </tr> <tr> <td>4</td> <td>91</td> <td>43</td> </tr> <tr> <td>5</td> <td>34</td> <td>16</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th>Hypertension</th> <th>Diabetes</th> <th>Renal failure</th> </tr> </thead> <tbody> <tr> <td>88%</td> <td>41%</td> <td>18%</td> </tr> <tr> <td>N = 187</td> <td>N = 87</td> <td>N = 38</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th colspan="3">Proximal anastomosis</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>Common femoral artery</td> <td>212</td> <td>100</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th colspan="3">Distal anastomosis</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>BK popliteal</td> <td>154</td> <td>73</td> </tr> <tr> <td>TP trunk</td> <td>58</td> <td>27</td> </tr> </tbody> </table> | Rutherford classification | | | | N | % | 3 | 87 | 41 | 4 | 91 | 43 | 5 | 34 | 16 | Hypertension | Diabetes | Renal failure | 88% | 41% | 18% | N = 187 | N = 87 | N = 38 | Proximal anastomosis | | | | N | % | Common femoral artery | 212 | 100 | Distal anastomosis | | | | N | % | BK popliteal | 154 | 73 | TP trunk | 58 | 27 | <ul style="list-style-type: none"> Retrospective, non-randomized, multicenter analysis Patients were considered at a low risk for thrombosis: Exclusion criteria included re-do operations, poor plantar and perimalleolar circulation, severe tissue loss and single-vessel runoff Secondary patency at 5 and 9 years was 58% and limb salvage at 5 and 9 years was 93% |
| Rutherford classification | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 87 | 41 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 91 | 43 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 34 | 16 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Hypertension | Diabetes | Renal failure | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 88% | 41% | 18% | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| N = 187 | N = 87 | N = 38 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Proximal anastomosis | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Common femoral artery | 212 | 100 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Distal anastomosis | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| BK popliteal | 154 | 73 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| TP trunk | 58 | 27 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

“...in low thrombotic risk patients mid-and long-term patency of vein and Propaten® graft [GORE® PROPATEN® Vascular Graft] is comparable. In case of PTFE use, we reported shorter surgery time, reduced hospital stay and wound complications. These observations led us to primarily choose the prosthetic graft in that subset of cases, saving the VSG [great saphenous vein] for distal revascularization in case of occlusive disease progression.” — V. Monaca


Results from an Italian multicentric registry comparing heparin-bonded ePTFE graft and autologous saphenous vein in below-knee femoro-popliteal bypasses⁵

Dorigo et al. 2012

Department of Vascular Surgery, University of Florence, Florence, Italy



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics | | Study details | | |
|--|---|--|--|----|
| Arterial hypertension 87% N = 482 | Diabetes 46% N = 253 | <ul style="list-style-type: none"> Retrospective, non-randomized, multicenter analysis Comparing GORE® PROPATEN® Vascular Graft and autologous saphenous vein There was no statistically significant difference in secondary patency or limb salvage between GORE® PROPATEN® Vascular Graft and autologous saphenous vein | | |
| Chronic renal failure 15% N = 82 | History of smoking 72% N = 403 | | | |
| | N | | | % |
| Hyperlipemia | 330 | | | 59 |
| Coronary artery disease | 251 | | | 45 |
| BK fem-pop | 414 | 75 | | |
| Infrapopliteal | 142 | 26 | | |
| TP trunk | 69 | 13 | | |
| Anterior tibial | 27 | 5 | | |
| Posterior tibial | 35 | 6 | | |
| Peroneal | 11 | 2 | | |

“...we had a 13% increase in secondary patency rates in ePTFE group, whereas the corresponding figure was only 6% in patients with occluded vein, thus confirming both the possibility of effectively treating occluded heparin-bonded grafts and the difficulty of dealing with occluded vein bypasses.” — *W. Dorigo*

“In patients with critical limb ischemia, the rates both of amputations at 4 years and of amputation-free survival were not different between autologous vein and heparin-bonded ePTFE... and this is an encouraging result, considering that limb salvage probably represents the main outcome in all these critical patients.” — *W. Dorigo*

Midterm results from a multicenter registry on the treatment of infrainguinal critical limb ischemia using a heparin-bonded ePTFE graft⁶

Pulli et al. 2010

Department of Vascular Surgery, University of Florence, Italy



Primary patency of GORE® PROPATEN® Vascular Graft

| | 1 year | 2 years | 3 years |
|--------------------------|--------|---------|---------|
| BK fem-pop N = 238 | 75% | 67% | 61% |
| Infrapopliteal N = 86 | 66% | 57% | 52% |

| Patient characteristics* | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---------------------------|----|--|--|---|---|---|-----|----|---|-----|----|---|----|----|---------------|--|--|--|---|---|-------|-----|----|-------|-----|----|--|---|---|-----------------------|-----|----|--------------------|-----|----|----------|-----|----|-----------------------|----|----|---------------------|--|--|-----------|----|----|----------|----|---|--------------------|---|---|------------------|----|---|--|
| <table border="1"> <thead> <tr> <th colspan="3">Rutherford classification</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>4</td> <td>230</td> <td>54</td> </tr> <tr> <td>5</td> <td>143</td> <td>34</td> </tr> <tr> <td>6</td> <td>52</td> <td>12</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th colspan="3">Vessel runoff</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0 – 1</td> <td>186</td> <td>44</td> </tr> <tr> <td>2 – 3</td> <td>239</td> <td>56</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>Arterial hypertension</td> <td>367</td> <td>86</td> </tr> <tr> <td>History of smoking</td> <td>326</td> <td>77</td> </tr> <tr> <td>Diabetes</td> <td>192</td> <td>45</td> </tr> <tr> <td>Chronic renal failure</td> <td>72</td> <td>17</td> </tr> <tr> <td colspan="3">Out-flow procedures</td> </tr> <tr> <td>Vein cuff</td> <td>45</td> <td>11</td> </tr> <tr> <td>Patching</td> <td>23</td> <td>5</td> </tr> <tr> <td>Tibial angioplasty</td> <td>5</td> <td>1</td> </tr> <tr> <td>Other procedures</td> <td>12</td> <td>3</td> </tr> </tbody> </table> | Rutherford classification | | | | N | % | 4 | 230 | 54 | 5 | 143 | 34 | 6 | 52 | 12 | Vessel runoff | | | | N | % | 0 – 1 | 186 | 44 | 2 – 3 | 239 | 56 | | N | % | Arterial hypertension | 367 | 86 | History of smoking | 326 | 77 | Diabetes | 192 | 45 | Chronic renal failure | 72 | 17 | Out-flow procedures | | | Vein cuff | 45 | 11 | Patching | 23 | 5 | Tibial angioplasty | 5 | 1 | Other procedures | 12 | 3 | <ul style="list-style-type: none"> Retrospective, non-randomized, multicenter study All patients had CLI Combined fem-pop and infrapopliteal limb salvage rate at 3 years was 81% in these CLI patients |
| Rutherford classification | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 230 | 54 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 143 | 34 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | 52 | 12 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Vessel runoff | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0 – 1 | 186 | 44 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 – 3 | 239 | 56 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Arterial hypertension | 367 | 86 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| History of smoking | 326 | 77 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Diabetes | 192 | 45 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Chronic renal failure | 72 | 17 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Out-flow procedures | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Vein cuff | 45 | 11 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Patching | 23 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tibial angioplasty | 5 | 1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Other procedures | 12 | 3 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

* Total N = 425.

“...the good results at 3 years in primary interventions in patients with more than one distal vessel and with rest pain could suggest a significant role of the heparin-bonded graft in these subgroups of patients.” — R. Pulli

“Primary and secondary patency rates make this graft an excellent alternative to autologous saphenous vein when it is absent, unsuitable, or of poor quality.” — R. Pulli

Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results⁷

Daenens et al. 2009
University Hospital Gasthuisberg, Belgium



Primary patency of GORE® PROPATEN® Vascular Graft

| | 1 year | 2 years |
|--------------------------|------------|------------|
| BK fem-pop N = 57 | 92% | 83% |
| Infrapopliteal N = 97 | 79% | 69% |

| Patient characteristics | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---------------|----|---|---|----|----|---|----|----|---|----|----|---|----|---|--|---|---|------------|----|----|----------------|----|----|----------|----|----|-----------------|----|----|------------------|----|----|----------|----|----|---|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>63</td> <td>26</td> </tr> <tr> <td>4</td> <td>60</td> <td>25</td> </tr> <tr> <td>5</td> <td>84</td> <td>35</td> </tr> <tr> <td>6</td> <td>17</td> <td>7</td> </tr> </tbody> </table> <p>Smokers 62% N = 149</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>BK fem-pop</td> <td>57</td> <td>37</td> </tr> <tr> <td>Infrapopliteal</td> <td>97</td> <td>63</td> </tr> <tr> <td>TP trunk</td> <td>10</td> <td>10</td> </tr> <tr> <td>Anterior tibial</td> <td>37</td> <td>38</td> </tr> <tr> <td>Posterior tibial</td> <td>29</td> <td>30</td> </tr> <tr> <td>Peroneal</td> <td>21</td> <td>22</td> </tr> </tbody> </table> | | N | % | 3 | 63 | 26 | 4 | 60 | 25 | 5 | 84 | 35 | 6 | 17 | 7 | | N | % | BK fem-pop | 57 | 37 | Infrapopliteal | 97 | 63 | TP trunk | 10 | 10 | Anterior tibial | 37 | 38 | Posterior tibial | 29 | 30 | Peroneal | 21 | 22 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center study Compared results from GORE® PROPATEN® Vascular Graft to autologous vein bypasses Adjunctive techniques: <ul style="list-style-type: none"> Below-knee fem-pop: 2 Miller cuff, 2 Taylor patch, 2 Linton patch Below-knee fem-distal: 11 Miller cuff, 3 Taylor patch, 7 Linton patch, 15 AV fistula The 2-year limb salvage rates for below-knee fem-pop and fem-distal bypasses were 98% and 87%, respectively 2-year autologous vein patency rates for below-knee fem-pop and fem-distal were 72% and 64%, respectively |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 63 | 26 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 60 | 25 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 84 | 35 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | 17 | 7 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| BK fem-pop | 57 | 37 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Infrapopliteal | 97 | 63 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| TP trunk | 10 | 10 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Anterior tibial | 37 | 38 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Posterior tibial | 29 | 30 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Peroneal | 21 | 22 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

“In this large retrospective study, heparin-bonded ePTFE grafts had 1- and 2-year primary patency results that were not significantly different from those for ASV grafts. Results in BK FP and FC applications were especially promising.” — *K. Daenens*

“Overall, our results...provide solid additional evidence that heparin-bonded ePTFE grafts represent an important new option in the treatment of peripheral arterial disease.”

— *K. Daenens*

PEPE II — A multicenter study with an end-point heparin-bonded expanded polytetrafluoroethylene vascular graft for above and below knee bypass surgery: determinants of patency⁸

Hugl et al. 2009

Department of Vascular Surgery, Medical University, Innsbruck, Austria



Primary patency of GORE® PROPATEN® Vascular Graft

1 year

| | |
|--------------------------|------------|
| BK fem-pop N = 37 | 74% |
| Infrapopliteal N = 15 | 79% |

| Patient characteristics* | | | Study details | | | | | | | | | | | | | | | |
|--|----------------------|----------------------|---------------|---|---|---------|----|-----|----------|----|----|-----------|----|----|----------|----|----|--|
| <p>Fontaine classification stage</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>Stage I</td> <td>1</td> <td>< 1</td> </tr> <tr> <td>Stage II</td> <td>62</td> <td>45</td> </tr> <tr> <td>Stage III</td> <td>26</td> <td>19</td> </tr> <tr> <td>Stage IV</td> <td>50</td> <td>36</td> </tr> </tbody> </table> | | | | N | % | Stage I | 1 | < 1 | Stage II | 62 | 45 | Stage III | 26 | 19 | Stage IV | 50 | 36 | <ul style="list-style-type: none"> • Prospective, non-randomized, multicenter study • Patients without suitable autologous vein • The 1-year secondary patency rates for below-knee fem-pop and infrapopliteal bypasses were 79% and 85%, respectively • Overall 1-year patency and limb salvage rates were 80% and 96%, respectively* |
| | N | % | | | | | | | | | | | | | | | | |
| Stage I | 1 | < 1 | | | | | | | | | | | | | | | | |
| Stage II | 62 | 45 | | | | | | | | | | | | | | | | |
| Stage III | 26 | 19 | | | | | | | | | | | | | | | | |
| Stage IV | 50 | 36 | | | | | | | | | | | | | | | | |
| <p>Vessel runoff</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>40</td> <td>29</td> </tr> <tr> <td>2</td> <td>50</td> <td>36</td> </tr> <tr> <td>3</td> <td>47</td> <td>34</td> </tr> </tbody> </table> | | | | N | % | 1 | 40 | 29 | 2 | 50 | 36 | 3 | 47 | 34 | | | | |
| | N | % | | | | | | | | | | | | | | | | |
| 1 | 40 | 29 | | | | | | | | | | | | | | | | |
| 2 | 50 | 36 | | | | | | | | | | | | | | | | |
| 3 | 47 | 34 | | | | | | | | | | | | | | | | |
| Hypertension | Diabetes | Current tobacco use | | | | | | | | | | | | | | | | |
| 71% N = 98 | 42% N = 59 | 42% N = 58 | | | | | | | | | | | | | | | | |

* Total N = 139, which includes 87 patients with above knee bypasses.

“...present data show that using the endpoint heparin bonded ePTFE graft for lower limb revascularization produces excellent results for AK bypasses and encouraging results for BK bypasses, when compared with data obtained from studies which used other prosthetic material. These encouraging results for BK bypasses were even seen in the subgroup of patients that generally have worse revascularization results due to the presence of risk factors.” — B. Hugl

“...our data suggests that the use of this graft is an excellent option when no autologous vein is available.” — B. Hugl

Heparin-bonded expanded polytetrafluoroethylene graft for infragenicular bypass: five-year results⁹

Lösel-Sadée & Alefelder. 2009

Department of Vascular Surgery, Sana Kliniken Dusseldorf, Dusseldorf, Germany



Primary patency of GORE® PROPATEN® Vascular Graft

| | 1 year | 2 years | 3 years | 4 years | 5 years |
|--------------------------|------------|------------|------------|------------|------------|
| BK fem-pop N = 30 | 77% | 71% | 71% | 71% | 71% |
| Infrapopliteal N = 45 | 64% | 57% | 50% | 50% | 50% |

| Patient characteristics* | | | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|----|-----|---------------|---|---|---|---|---|---|----|----|---|----|----|---|---|-----|--|---|---|---|----|----|---|----|----|---|----|----|--|---|---|-------------|----|----|----------------|----|----|----------|----|----|-----------------|----|----|------------------|---|---|----------|---|---|---------------|---|-----|--|--|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>6</td> <td>8</td> </tr> <tr> <td>4</td> <td>43</td> <td>57</td> </tr> <tr> <td>5</td> <td>25</td> <td>33</td> </tr> <tr> <td>6</td> <td>1</td> <td>< 1</td> </tr> </tbody> </table> <p>Vessel runoff</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>41</td> <td>55</td> </tr> <tr> <td>2</td> <td>18</td> <td>24</td> </tr> <tr> <td>3</td> <td>16</td> <td>21</td> </tr> </tbody> </table> <p>Hypertension 96% (N = 72) Diabetes 49% (N = 37) Renal insufficiency 31% (N = 23)</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>BK fem-pop*</td> <td>30</td> <td>40</td> </tr> <tr> <td>Infrapopliteal</td> <td>45</td> <td>60</td> </tr> <tr> <td>TP trunk</td> <td>12</td> <td>27</td> </tr> <tr> <td>Anterior tibial</td> <td>14</td> <td>31</td> </tr> <tr> <td>Posterior tibial</td> <td>9</td> <td>2</td> </tr> <tr> <td>Peroneal</td> <td>9</td> <td>2</td> </tr> <tr> <td>Dosalis pedis</td> <td>1</td> <td>< 1</td> </tr> </tbody> </table> | | | | N | % | 3 | 6 | 8 | 4 | 43 | 57 | 5 | 25 | 33 | 6 | 1 | < 1 | | N | % | 1 | 41 | 55 | 2 | 18 | 24 | 3 | 16 | 21 | | N | % | BK fem-pop* | 30 | 40 | Infrapopliteal | 45 | 60 | TP trunk | 12 | 27 | Anterior tibial | 14 | 31 | Posterior tibial | 9 | 2 | Peroneal | 9 | 2 | Dosalis pedis | 1 | < 1 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center study Vein cuffs were created at the distal anastomosis in 5 patients; no patches were used The 3- and 4-year secondary patency rates for below-knee fem-pop and infrapopliteal bypasses were 83% and 72%, respectively The 5-year limb salvage rate was 84% <p>“The primary patency results are especially encouraging in light of the fact that the patients were seriously ill, as indicated by the high rates of Rutherford category 4 to 6 disease, renal insufficiency, and previous treatment for PAD and the low rate of multiple-vessel runoff in the series.” — H. Lösel-Sadée</p> | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 6 | 8 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 43 | 57 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 25 | 33 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | 1 | < 1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | 41 | 55 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | 18 | 24 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 16 | 21 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| BK fem-pop* | 30 | 40 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Infrapopliteal | 45 | 60 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| TP trunk | 12 | 27 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Anterior tibial | 14 | 31 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Posterior tibial | 9 | 2 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Peroneal | 9 | 2 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Dosalis pedis | 1 | < 1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

* Distal of the knee articulation.

Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass?¹⁰

Peeters et al. 2008

Department of Cardiovascular and Thoracic Surgery, Imelda Hospital, Bonheiden, Belgium



Primary patency of GORE® PROPATEN® Vascular Graft


| | 1 year | 2 years | 3 years |
|--------------------------|------------|------------|------------|
| BK fem-pop N = 41 | 86% | 79% | 75% |
| Infrapopliteal N = 37 | 71% | 60% | 60% |

| Patient characteristics* | | | Study details | | | | | | | | | | | | | | | | | | | | |
|--|---|---|---------------|---|---|------------|----|----|----------------|----|----|----------|----|----|---|----|----|------------------|----|----|----------|---|----|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>85</td> <td>56</td> </tr> <tr> <td>4</td> <td>29</td> <td>19</td> </tr> <tr> <td>5</td> <td>39</td> <td>26</td> </tr> </tbody> </table> | | | | N | % | 3 | 85 | 56 | 4 | 29 | 19 | 5 | 39 | 26 | <ul style="list-style-type: none"> • Prospective, non-randomized, multicenter study • No adjunctive techniques (patches or cuffs) were used in the study • The 3-year secondary patency rate for below-knee fem-pop and infrapopliteal bypasses were 80% and 62%, respectively • The 3-year limb salvage rate for all CLI patients was 86%* | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | |
| 3 | 85 | 56 | | | | | | | | | | | | | | | | | | | | | |
| 4 | 29 | 19 | | | | | | | | | | | | | | | | | | | | | |
| 5 | 39 | 26 | | | | | | | | | | | | | | | | | | | | | |
| <p>Vessel runoff</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>7</td> <td>5</td> </tr> <tr> <td>1</td> <td>76</td> <td>50</td> </tr> <tr> <td>2</td> <td>43</td> <td>28</td> </tr> <tr> <td>3</td> <td>27</td> <td>18</td> </tr> </tbody> </table> | | | | N | % | 0 | 7 | 5 | 1 | 76 | 50 | 2 | 43 | 28 | 3 | 27 | 18 | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | |
| 0 | 7 | 5 | | | | | | | | | | | | | | | | | | | | | |
| 1 | 76 | 50 | | | | | | | | | | | | | | | | | | | | | |
| 2 | 43 | 28 | | | | | | | | | | | | | | | | | | | | | |
| 3 | 27 | 18 | | | | | | | | | | | | | | | | | | | | | |
| <p>Hypertension</p> <p>64%</p> <p>N = 88</p> | <p>Diabetes</p> <p>27%</p> <p>N = 37</p> | <p>Renal insufficiency</p> <p>9%</p> <p>N = 12</p> | | | | | | | | | | | | | | | | | | | | | |
| <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>BK fem-pop</td> <td>41</td> <td>53</td> </tr> <tr> <td>Infrapopliteal</td> <td>37</td> <td>47</td> </tr> <tr> <td>TP trunk</td> <td>11</td> <td>30</td> </tr> <tr> <td>Anterior tibial</td> <td>10</td> <td>27</td> </tr> <tr> <td>Posterior tibial</td> <td>10</td> <td>27</td> </tr> <tr> <td>Peroneal</td> <td>6</td> <td>16</td> </tr> </tbody> </table> | | | | N | % | BK fem-pop | 41 | 53 | Infrapopliteal | 37 | 47 | TP trunk | 11 | 30 | Anterior tibial | 10 | 27 | Posterior tibial | 10 | 27 | Peroneal | 6 | 16 |
| | N | % | | | | | | | | | | | | | | | | | | | | | |
| BK fem-pop | 41 | 53 | | | | | | | | | | | | | | | | | | | | | |
| Infrapopliteal | 37 | 47 | | | | | | | | | | | | | | | | | | | | | |
| TP trunk | 11 | 30 | | | | | | | | | | | | | | | | | | | | | |
| Anterior tibial | 10 | 27 | | | | | | | | | | | | | | | | | | | | | |
| Posterior tibial | 10 | 27 | | | | | | | | | | | | | | | | | | | | | |
| Peroneal | 6 | 16 | | | | | | | | | | | | | | | | | | | | | |
| <p>“Propaten® Vascular Graft [GORE® PROPATEN® Vascular Graft] may succeed in bridging the gap between venous and regular PTFE bypass...especially for infrapopliteal bypasses.”</p> <p>— P. Peeters</p> | | | | | | | | | | | | | | | | | | | | | | | |

* Limb salvage rates are for both above knee and below-knee bypasses.

Overall weighted average* primary patency in above-knee bypasses†

| | 1 year | 2 years | 3 years | 4 years | 5 years |
|---|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Vein | 76% N = 242 | 70% N = 236 | 71% N = 109 | ‡ | 60% N = 95 |
| GORE® PROPATEN® Vascular Graft | 84% N = 860 | 79% N = 606 | 76% N = 520 | 74% N = 445 | 68% N = 445 |



* Weighted Average =
$$\frac{(N_1 \times \text{Primary Patency}_1) + (N_2 \times \text{PP}_2) + \dots + (N_n \times \text{PP}_n)}{N_1 + N_2 + \dots + N_n}$$

† Above-knee (AK) inclusion criteria for GORE® PROPATEN® Vascular Graft literature used in this analysis were (1) articles in English language, (2) clinical journal articles or book chapters, (3) non-overlapping patient populations and (4) AK bypass primary patency reported for at least 12 months of follow-up. Additional exclusion criteria were (1) reviews, case reports or meta-analysis articles and (2) articles containing the key word AV access (including synonyms). (data on file 2019; W. L. Gore & Associates, Inc; Flagstaff, AZ.)

‡ No data available.


Results from a multicenter registry of heparin-bonded expanded polytetrafluoroethylene graft for above-the-knee femoropopliteal bypass?¹¹

Piffaretti et al. 2018

Università degli studi dell'Insubria, Varese, Italy



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics | Study details | | | | | | | | | | | | | | | |
|--|---------------|----|---|---|-----|----|---|----|----|---|----|----|---|----|---|--|
| <p>Rutherford classification</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>200</td> <td>55</td> </tr> <tr> <td>4</td> <td>86</td> <td>24</td> </tr> <tr> <td>5</td> <td>68</td> <td>18</td> </tr> <tr> <td>6</td> <td>10</td> <td>3</td> </tr> </tbody> </table> <p>Hypertension 81% N = 294</p> <p>Diabetes 39% N = 141</p> <p>History of smoking 59% N = 216</p> | | N | % | 3 | 200 | 55 | 4 | 86 | 24 | 5 | 68 | 18 | 6 | 10 | 3 | <ul style="list-style-type: none"> Retrospective, non-randomized, multicenter analysis At the time of publication, the GORE® PROPATEN® Vascular Graft Italian registry had tracked 1,401 interventions performed for peripheral arterial obstructive disease using GORE® PROPATEN® Vascular Graft in a “real-world” setting Estimated survival at 5 years was 75% Estimated amputation-free survival at 5 years was 74% Estimated assisted primary patency, secondary patency and limb salvage at 5 years were 65%, 75% and 95%, respectively Postoperative medical treatment with warfarin alone was found to be an independent risk factor for loss of primary patency compared to dual antiplatelet therapy |
| | N | % | | | | | | | | | | | | | | |
| 3 | 200 | 55 | | | | | | | | | | | | | | |
| 4 | 86 | 24 | | | | | | | | | | | | | | |
| 5 | 68 | 18 | | | | | | | | | | | | | | |
| 6 | 10 | 3 | | | | | | | | | | | | | | |

“[Above-the-knee femoropopliteal bypass] with the use of HB-ePTFE remains an effective option, with low rate of perioperative complications and satisfactory long-term results.”

— G. Piffaretti

“...in our opinion, “[Above-the-knee femoropopliteal bypass] is a valid and viable first-line alternative to endovascular surgery in long or complex lesions of the SFA.” — G. Piffaretti


Heparin-bonded expanded polytetrafluoroethylene femoropopliteal bypass grafts outperform expanded polytetrafluoroethylene grafts without heparin in a long-term comparison¹²

Samson et al. 2016

Sarasota Vascular Specialists in Sarasota, Florida, USA



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | | | Study details |
|--------------------------|--------|--------|--|
| | N | % | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center review of prospectively collected data No adjunctive technique (patches or cuffs) were used Most grafts were 6 mm ring reinforced PLAVIX® Clopidogrel Bisulfate usage had a significant benefit on overall primary patency Loss of patency was found to be related to younger age, absence of claudication, isolated popliteal artery and lower post-operative ABI As early as 3 months, a significant difference in patency favoring GORE® PROPATEN® Vascular Graft was seen and was maintained at 5 years (75% versus 56%) and in both AK (85% versus 59%) and BK (60% versus 0%/undeterminable) locations |
| Gangrene | 27 | 20 | |
| Ulceration | 39 | 30 | |
| Rest pain | 40 | 30 | |
| Claudication | 43 | 33 | |
| Runoff vessels* | | | |
| | N | % | |
| 0 | 2 | 2 | |
| 1 | 44 | 37 | |
| 2 | 49 | 41 | |
| 3 | 24 | 20 | |
| Hypertension | 84% | | |
| Diabetes | 45% | | |
| Smoking | 32% | | |
| N = 110 | N = 59 | N = 42 | |

* GORE® PROPATEN® Vascular Graft group.

“These data show that the Propaten HePTFE graft [GORE® PROPATEN® Vascular Graft] offered significantly better long-term patency over the SePTFE graft, suggesting Propaten [GORE® PROPATEN® Vascular Graft] as the prosthetic graft of choice for bypasses to the femoropopliteal artery when autologous vein is unavailable or inappropriate.” — R. Samson

“Propaten HePTFE grafts [GORE® PROPATEN® Vascular Grafts] seem to be as effective as vein for AK femoropopliteal artery bypass. Because of the excellent results observed in this series, we now use the Propaten graft [GORE® PROPATEN® Vascular Graft] preferentially over great saphenous vein for AK bypass except in younger patients with available appropriate autologous conduit.” — R. Samson

The Scandinavian Propaten® Trial — 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses — a randomized clinical controlled multi-centre trial¹³

Lindholt et al. 2011

Vascular Research Unit, Department of Vascular Surgery, Viborg Hospital, Denmark



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | Study details |
|--|--|
| <p>Smokers</p> <p>53%</p> <p>N = 144</p> | <ul style="list-style-type: none"> • Prospective, randomized, multicenter (11 centers) study comparing GORE® PROPATEN® Vascular Graft versus standard ePTFE • Fem-pop (majority above-knee) and fem-fem bypasses • Statistically significant improvement in primary and secondary patency with GORE® PROPATEN® Vascular Graft versus standard ePTFE for all bypasses • In fem-pop patients with CLI, GORE® PROPATEN® Vascular Graft primary patency was 80% while standard ePTFE patency was 58% ($P < 0.05$) • GORE® PROPATEN® Vascular Graft reduced the risk of graft occlusion by 40% overall and by 50% in patients with CLI |
| <p>Diabetes</p> <p>15%</p> <p>N = 39</p> | |
| <p>Critical limb ischemia (CLI)</p> <p>36%</p> <p>N = 100</p> | |
| <p>* GORE® PROPATEN® Vascular Graft group, N = 272, which includes 160 patients with fem-fem bypass.</p> | |

“We have seen that the GORE® PROPATEN® Vascular Graft keeps its promise as shown in previously conducted prospective and retrospective studies.”

— J. Lindholt

Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results⁷

Daenens et al. 2009

University Hospital Gasthuisberg, Belgium



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | Study details | | | | | | | | | | | | | | | | | | |
|--|---------------------------|----|--|--|---|---|---|----|----|---|----|----|---|----|----|---|----|---|--|
| <table border="1"> <thead> <tr> <th colspan="3">Rutherford classification</th> </tr> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>63</td> <td>26</td> </tr> <tr> <td>4</td> <td>60</td> <td>25</td> </tr> <tr> <td>5</td> <td>84</td> <td>35</td> </tr> <tr> <td>6</td> <td>17</td> <td>7</td> </tr> </tbody> </table> <p>Smokers 62% N = 149</p> <p>Redo bypass 36% N = 86</p> | Rutherford classification | | | | N | % | 3 | 63 | 26 | 4 | 60 | 25 | 5 | 84 | 35 | 6 | 17 | 7 | <ul style="list-style-type: none"> Retrospective, non-randomized, single-center study Study compared results from GORE® PROPATEN® Vascular Graft to autologous vein bypasses The 1- and 2-year primary patency rates for above-knee fem-pop bypasses using autologous vein were 91% and 80%, respectively |
| Rutherford classification | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | |
| 3 | 63 | 26 | | | | | | | | | | | | | | | | | |
| 4 | 60 | 25 | | | | | | | | | | | | | | | | | |
| 5 | 84 | 35 | | | | | | | | | | | | | | | | | |
| 6 | 17 | 7 | | | | | | | | | | | | | | | | | |

* GORE® PROPATEN® Vascular Graft group, N = 240.

“Overall, our results...provide solid additional evidence that heparin-bonded ePTFE grafts represent an important new option in the treatment of peripheral arterial disease.”

— K. Daenens


Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass?¹⁰

Peeters et al. 2008

Department of Cardiovascular and Thoracic Surgery, Imelda Hospital, Bonheiden, Belgium



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|--|----|---|---|----|----|---|----|----|---|----|----|--|---|---|---|---|---|---|----|----|---|----|----|---|----|----|--|--|--|--|---|
| <p>Rutherford classification*</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>3</td> <td>85</td> <td>56</td> </tr> <tr> <td>4</td> <td>29</td> <td>19</td> </tr> <tr> <td>5</td> <td>39</td> <td>25</td> </tr> </tbody> </table> <p>Runoff vessels*</p> <table border="1"> <thead> <tr> <th></th> <th>N</th> <th>%</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>7</td> <td>5</td> </tr> <tr> <td>1</td> <td>76</td> <td>50</td> </tr> <tr> <td>2</td> <td>43</td> <td>28</td> </tr> <tr> <td>3</td> <td>27</td> <td>18</td> </tr> </tbody> </table> <table border="0"> <tr> <td> <p>Diabetes*</p> <p>27%</p> <p>N = 37</p> </td> <td> <p>Hypertension†</p> <p>64%</p> <p>N = 88</p> </td> </tr> <tr> <td> <p>Coronary artery disease</p> <p>30%</p> <p>N = 42</p> </td> <td> <p>Nicotine use†</p> <p>48%</p> <p>N = 66</p> </td> </tr> </table> | | N | % | 3 | 85 | 56 | 4 | 29 | 19 | 5 | 39 | 25 | | N | % | 0 | 7 | 5 | 1 | 76 | 50 | 2 | 43 | 28 | 3 | 27 | 18 | <p>Diabetes*</p> <p>27%</p> <p>N = 37</p> | <p>Hypertension†</p> <p>64%</p> <p>N = 88</p> | <p>Coronary artery disease</p> <p>30%</p> <p>N = 42</p> | <p>Nicotine use†</p> <p>48%</p> <p>N = 66</p> | <ul style="list-style-type: none"> • Prospective, non-randomized, multicenter study • No adjunctive techniques (patches or cuffs) were used in the study • 73% of patients had been previously treated for peripheral vascular disease (PVD) • The 3-year limb salvage rate for all CLI patients was 86%‡ |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 85 | 56 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 29 | 19 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 39 | 25 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0 | 7 | 5 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | 76 | 50 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | 43 | 28 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 27 | 18 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>Diabetes*</p> <p>27%</p> <p>N = 37</p> | <p>Hypertension†</p> <p>64%</p> <p>N = 88</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <p>Coronary artery disease</p> <p>30%</p> <p>N = 42</p> | <p>Nicotine use†</p> <p>48%</p> <p>N = 66</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

* Total N = 153 limbs.

† Total N = 138, which includes 97 patients that underwent below-knee bypasses.

‡ Limb salvage rates are for both above and below-knee bypasses.

“...it is our opinion that the Propaten Vascular Graft [GORE® PROPATEN® Vascular Graft] may succeed in bridging the gap between venous conduits and regular ePTFE grafts.”

— P. Peeters


Heparin-bonded expanded polytetrafluoroethylene vascular graft for femoropopliteal and femorocrural bypass grafting: 1-year results¹⁴

Bosiers et al. 2006

Dendermonde, Bonheiden, and Genk, Belgium



Primary patency of GORE® PROPATEN® Vascular Graft



| Patient characteristics* | Study details | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|---------------|----------------------|---|---|----|----|---|----|----|---|----|----|--|---|---|---|---|---|---|----|----|---|----|----|---|----|----|---------------|-----------|----------------------|------------|------------|-----------|--------|--------|-------|---|
| <p style="text-align: center;">Rutherford classification*</p> <table border="1" style="margin-left: auto; margin-right: auto; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">N</th> <th style="text-align: center;">%</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">58</td> <td style="text-align: center;">58</td> </tr> <tr> <td style="text-align: center;">4</td> <td style="text-align: center;">16</td> <td style="text-align: center;">16</td> </tr> <tr> <td style="text-align: center;">5</td> <td style="text-align: center;">26</td> <td style="text-align: center;">26</td> </tr> </tbody> </table> <p style="text-align: center;">Runoff vessels*</p> <table border="1" style="margin-left: auto; margin-right: auto; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">N</th> <th style="text-align: center;">%</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">0</td> <td style="text-align: center;">7</td> <td style="text-align: center;">7</td> </tr> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">51</td> <td style="text-align: center;">51</td> </tr> <tr> <td style="text-align: center;">2</td> <td style="text-align: center;">21</td> <td style="text-align: center;">21</td> </tr> <tr> <td style="text-align: center;">3</td> <td style="text-align: center;">21</td> <td style="text-align: center;">21</td> </tr> </tbody> </table> <table style="width: 100%; text-align: center;"> <tr> <td style="width: 33%;">Hypertension†</td> <td style="width: 33%;">Diabetes†</td> <td style="width: 33%;">Renal insufficiency†</td> </tr> <tr> <td style="font-size: 2em; color: red;">72%</td> <td style="font-size: 2em; color: red;">28%</td> <td style="font-size: 2em; color: red;">9%</td> </tr> <tr> <td>N = 62</td> <td>N = 24</td> <td>N = 8</td> </tr> </table> | | N | % | 3 | 58 | 58 | 4 | 16 | 16 | 5 | 26 | 26 | | N | % | 0 | 7 | 7 | 1 | 51 | 51 | 2 | 21 | 21 | 3 | 21 | 21 | Hypertension† | Diabetes† | Renal insufficiency† | 72% | 28% | 9% | N = 62 | N = 24 | N = 8 | <ul style="list-style-type: none"> Prospective, non-randomized, multicenter study All grafts were thin-wall, ringed, 6 mm diameter configurations Angiography performed immediately after graft implantation to detect technical failures Postoperatively, patients received PLAVIX® Clopidogrel Bisulfate usage (75 mg/day) for the first month after surgery, low-molecular-weight heparin (0.6 mL/day) for the first three weeks and ASPIRIN Acetylsalicylic Acid (100 mg/day) indefinitely 96% secondary patency at one-year for above-knee fem-pop bypass |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 58 | 58 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | 16 | 16 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | 26 | 26 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | N | % | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0 | 7 | 7 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | 51 | 51 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | 21 | 21 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | 21 | 21 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Hypertension† | Diabetes† | Renal insufficiency† | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 72% | 28% | 9% | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| N = 62 | N = 24 | N = 8 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

* N = 100 limbs.

† N = 86 patients.


“In the light of these in vivo results, we speculate that a decrease in platelet and thrombus deposition on the CBAS ePTFE graft surface may have contributed to the promising 1-year patency rates in our clinical series. Although amelioration of intimal hyperplasia is not the primary target of heparinization technology, it is intriguing to consider the possibility that the CBAS graft surface might simultaneously address two sources of graft failure: thrombosis and intimal hyperplasia.” — M. Bosiers

“Our findings... indicate that use of this graft is an excellent option for infrainguinal bypass grafting in patients with peripheral vascular disease when autologous vein is not available.” — M. Bosiers

References

1. Kaisar J, Chen A, Cheung M, Kfoury E, Bechara C, Lin P. Comparison of propaten heparin-bonded vascular graft with distal anastomotic patch versus autogenous saphenous vein graft in tibial artery bypass. *Vascular* 2018;26(2):117-125.
2. Uhl C, Hock C, Betz T, Greindl M, Töpel I, Steinbauer M. Comparison of venous and HePTFE tibial and peroneal bypasses in critical limb ischemia patients unsuitable for endovascular revascularization. *Vascular* 2015;23(6):607-613.
3. Neville RF, Babrowicz J, Amdur R, Sidawy A. Heparin-bonded ePTFE (Propaten): is it as good as autologous vein for tibial bypass? *Italian Journal of Vascular & Endovascular Surgery* 2014;21(1)49-56.
4. Monaca V, Battaglia G, Turiano SA, Tringale R, Catalfamo S. Subpopliteal revascularization. Criteria analysis for use of E-PTFE (Propaten®) as first choice conduit. *Italian Journal of Vascular & Endovascular Surgery* 2013;20(3):165-169.
5. Dorigo W, Raffaele P, Piffaretti G. Results from an Italian multicentric registry comparing heparin-bonded ePTFE graft and autologous saphenous vein in below-knee femoro-popliteal bypasses. *Journal of Cardiovascular Surgery* 2012;53(2):187-193.
6. Pulli R, Dorigo W, Castelli P, et al. Propaten Italian Registry Group. Midterm results from a multicenter registry on the treatment of infrainguinal critical limb ischemia using a heparin-bonded ePTFE graft. *Journal of Vascular Surgery* 2010;51(5):1167-1177.e1.
7. Daenens K, Schepers S, Fourneau I, Houthoofd S, Nevelsteen A. Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results. *Journal of Vascular Surgery* 2009;49(5):1210-1216.
8. Hugl B, Nevelsteen A, Daenens K, et al. PEPE II Study Group. PEPE II - a multicenter study with an end-point heparin-bonded expanded polytetrafluoroethylene vascular graft for above and below knee bypass surgery: determinants of patency. *Journal of Cardiovascular Surgery* 2009;50(2):195-203.
9. Lösel-Sadée H, Alefelder C. Heparin-bonded expanded polytetrafluoroethylene graft for infragenicular bypass: five-year results. *Journal of Cardiovascular Surgery* 2009;50(3):339-343.
10. Peeters P, Verbist J, Deloose K, Bosiers M. Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass? *Italian Journal of Vascular & Endovascular Surgery* 2008;15(3):143-148.
11. Piffaretti G, Dorigo W, Castelli P, Pratesi C, Pulli R; PROPATEN Italian Registry Group. Results from a multicenter registry of heparin-bonded expanded polytetrafluoroethylene graft for above-the-knee femoropopliteal bypass. *Journal of Vascular Surgery* 2018;67(5):1463-1471.e1.
12. Samson RH, Morales R, Showalter DP, Lepore MR, Nair DG. Heparin-bonded expanded polytetrafluoroethylene femoropopliteal bypass grafts outperform expanded polytetrafluoroethylene grafts without heparin in a long-term comparison. *Journal of Vascular Surgery* 2016;64(3):638-647.
13. Lindholt JS, Gottschalksen B, Johannesen N, et al. The Scandinavian Propaten® trial - 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses - a randomised clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673. <https://www.sciencedirect.com/science/article/pii/S1078588411000621>
14. Bosiers M, Deloose K, Verbist J, et al. Heparin-bonded expanded polytetrafluoroethylene vascular graft for femoropopliteal and femorocrural bypass grafting: 1-year results. *Journal of Vascular Surgery* 2006;43(2):313-319.
15. CBAS Heparin Surface. W. L. Gore & Associates website. Accessed April 14, 2020. <https://www.goremedical.com/cbas/references>.

Features Gore's CBAS Heparin Surface, the proven heparin bonding technology for lasting thromboresistance, used in many of Gore's interventional and vascular surgery products. End-point covalent bonding keeps heparin anchored to the graft surface, while the bioactive site remains free to interact with the blood to help prevent clotting.¹⁵



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific)
1800 680 424 (Australia/New Zealand)
00800 6334 4673 (Europe)

800 437 8181 (United States)
928 779 2771 (United States)

goremedical.com

 Consult Instructions
for Use
eifu.goremedical.com

INDICATIONS FOR USE: GORE® PROPATEN® Vascular Grafts are intended for use as vascular prostheses for replacement or bypass of diseased vessels in patients suffering occlusive or aneurysmal diseases, in trauma patients requiring vascular replacement, for dialysis access, or for other vascular procedures. **CONTRAINDICATIONS:** A. DO NOT use the GORE® PROPATEN® Vascular Graft in patients with known hypersensitivity to heparin, including those patients who have had a previous incidence of Heparin-Induced Thrombocytopenia (HIT) type II. B. DO NOT use any configuration of GORE® PROPATEN® Vascular Grafts with Removable Rings, Non-Removable Rings or Integrated Rings for coronary artery bypass or cerebral reconstruction procedures. C. DO NOT use GORE® PROPATEN® Vascular Grafts as a patch. If cut and used as a patch, GORE® PROPATEN® Vascular Grafts may lack adequate transverse strength. Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. ~~R~~only The articles listed here may include outcome descriptions that reflect use of devices outside the approved *Instructions for Use* (IFU), including data registries intended to provide broader information about patient outcomes. Gore does not promote the off-label use of devices. Physicians should consult the IFU for complete device information, including contraindications, warnings and cautions.

Products listed may not be available in all markets.

BAYER and ASPIRIN are trademarks of Bayer Intellectual Property GmbH or its affiliated companies.
MEDLINE is a trademark of Medline Industries, Inc.

SANOFI-AVENTIS and PLAVIX are trademarks of Sanofi-Aventis.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

GORE, PROPATEN and designs are trademarks of W. L. Gore & Associates.

© 2014, 2019, 2020 W. L. Gore & Associates, Inc. AT0145-EN5 APRIL 2020