



GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

CONFIGURATIONS



Catalogue number	Stent labeled/nominal diameter (mm)	Max. post-dilated stent diameter (mm)*	Crimped stent length (mm)	Catheter length (cm)	Guidewire diameter	Introducer sheath size (Fr)
BXA051502W	5	8	15	135	.035"	7
BXA051902W	5	8	19	135	.035"	7
BXA052902W	5	8	29	135	.035"	7
BXA053902W	5	8	39	135	.035"	7
BXA055902W	5	8	59	135	.035"	7
BXA057902W	5	8	79	135	.035"	7
BXA061502W	6	8	15	135	.035"	7
BXA061902W	6	8	19	135	.035"	7
BXA062902W	6	8	29	135	.035"	7
BXA063902W	6	8	39	135	.035"	7
BXA065902W	6	8	59	135	.035"	7
BXA067902W	6	8	79	135	.035"	7
BXA071502W	7	11	15	135	.035"	7
BXA071902W	7	11	19	135	.035"	7
BXA072902W	7	11	29	135	.035"	7
BXA073902W	7	11	39	135	.035"	7
BXA075902W	7	11	59	135	.035"	7
BXA077902W	7	11	79	135	.035"	7
BXA082902W	8	11	29	135	.035"	7
BXA083902W	8	11	39	135	.035"	7
BXA085902W	8	11	59	135	.035"	7
BXA087902W	8	11	79	135	.035"	8
BXA092902W	9	13	29	135	.035"	8
BXA093902W	9	13	39	135	.035"	8
BXA095902W	9	13	59	135	.035"	8
BXA097902W	9	13	79	135	.035"	8
BXA102902W	10	13	29	135	.035"	8
BXA103902W	10	13	39	135	.035"	8
BXA105902W	10	13	59	135	.035"	8
BXA107902W	10	13	79	135	.035"	8
BXA112902W	11	16 [†]	29	135	.035"	8
BXA113902W	11	16 [†]	39	135	.035"	8
BXA115902W	11	16 [†]	59	135	.035"	8
BXA117902W	11	16 [†]	79	135	.035"	8
BXAL082902W	8	16 [†]	29	135	.035"	7
BXAL083902W	8	16 [†]	39	135	.035"	7
BXAL085902W	8	16 [†]	59	135	.035"	8
BXAL087902W	8	16 [†]	79	135	.035"	8

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

† Technical limit of the device as determined by in-vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the Gore VBX FLEX Clinical Study.

Refer to *Instructions for Use* for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. ^{FX} only

Products listed may not be available in all markets.

GORE, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2021 W. L. Gore & Associates, Inc. 21393952-EN DECEMBER 2021



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65 67332882 (Asia Pacific)

1800 680 424 (Australia/New Zealand)

800 437 8181 (United States)

goremedical.com

0800 6334 4673 (Europe)

928 779 2771 (United States)