

SIZING CARD

GORE® VIABAHN® VBX

Balloon Expandable Endoprosthesis



Together, improving life



Device sizing

Introducer sheath (Fr)

Stent labeled/ nominal diameter (mm)	Crimped stent length						Maximum post-dilated stent diameter (mm)*
	15 mm	19 mm	29 mm	39 mm	59 mm	79 mm	
5	7	7	7	7	7	7	8
6	7	7	7	7	7	7	8
7	7	7	7	7	7	7	11
8	—	—	7	7	7	8	11
8L	—	—	7	7	8	8	16 [†]
9	—	—	8	8	8	8	13
10	—	—	8	8	8	8	13
11	—	—	8	8	8	8	16 [†]

All configurations available in .035" over-the-wire 135 cm catheter length.

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

† Technical limit of the device. Expansion beyond 13 mm is outside of the indication — see *Instructions for Use*.

Refer to *Instructions for Use* for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2022 W. L. Gore & Associates, Inc. 21407150-EN JANUARY 2022

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882
Australia/New Zealand 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

