

Device Selection



GORE® VIATORR® TIPS Endoprosthesis with Controlled Expansion




CATALOGUE NUMBER	ENDOPROSTHESIS INTERNAL DIAMETER (mm)*	GRAFT LINED LENGTH (cm)	GRAFT UNLINED LENGTH (cm)	MAXIMUM GUIDEWIRE DIAMETER (in)	HEMOSTATIC INTRODUCER SHEATH (Fr)	MAXIMUM DILATATION BALLOON DIAMETER (mm)*
PT8104275	8–10	4	2	≤ 0.035	10	10
PT8105275	8–10	5	2	≤ 0.035	10	10
PT8106275	8–10	6	2	≤ 0.035	10	10
PT8107275	8–10	7	2	≤ 0.035	10	10
PT8108275	8–10	8	2	≤ 0.035	10	10

* The selected balloon diameter may be 8, 9, or 10 mm, and should not exceed 10 mm. A balloon which can reach an inflation pressure of 10 ATM must be selected, and should be inflated to a minimum of 10 ATM.



PERFORMANCE
by design

Device Selection



GORE® VIATORR® TIPS Endoprosthesis

CATALOGUE NUMBER	ENDOPROSTHESIS INTERNAL DIAMETER (mm)	GRAFT LINED LENGTH (cm)	GRAFT UNLINED LENGTH (cm)	MAXIMUM GUIDEWIRE DIAMETER (in)	HEMOSTATIC INTRODUCER SHEATH (Fr)	MAXIMUM DILATION BALLOON DIAMETER (mm)
PT084275	8	4	2	≤ 0.038	10	8
PT085275	8	5	2	≤ 0.038	10	8
PT086275	8	6	2	≤ 0.038	10	8
PT087275	8	7	2	≤ 0.038	10	8
PT088275	8	8	2	≤ 0.038	10	8
PT104275	10	4	2	≤ 0.038	10	10
PT105275	10	5	2	≤ 0.038	10	10
PT106275	10	6	2	≤ 0.038	10	10
PT107275	10	7	2	≤ 0.038	10	10
PT108275	10	8	2	≤ 0.038	10	10

W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.

GORE®, PERFORMANCE BY DESIGN, VIATORR®, and designs are trademarks of W. L. Gore & Associates.
© 2016 W. L. Gore & Associates, Inc. AV0284-EN1 APRIL 2016



PERFORMANCE
by design