



GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

CONFIGURATIONS



Catalogue number	Stent labeled / nominal diameter (mm)	Max. post-dilated stent diameter (mm)*	Crimped stent length (mm)	Catheter length (cm)	Guidewire diameter	Introducer sheath size (Fr)
BXA051502B	5	8	15	135	0.035"	7
BXA051902B	5	8	19	135	0.035"	7
BXA052902B	5	8	29	135	0.035"	7
BXA053902B	5	8	39	135	0.035"	7
BXA055902B	5	8	59	135	0.035"	7
BXA057902B	5	8	79	135	0.035"	7
BXA061502B	6	8	15	135	0.035"	7
BXA061902B	6	8	19	135	0.035"	7
BXA062902B	6	8	29	135	0.035"	7
BXA063902B	6	8	39	135	0.035"	7
BXA065902B	6	8	59	135	0.035"	7
BXA067902B	6	8	79	135	0.035"	7
BXA071502B	7	11	15	135	0.035"	7
BXA071902B	7	11	19	135	0.035"	7
BXA072902B	7	11	29	135	0.035"	7
BXA073902B	7	11	39	135	0.035"	7
BXA075902B	7	11	59	135	0.035"	7
BXA077902B	7	11	79	135	0.035"	7
BXA082902B	8	11	29	135	0.035"	7
BXA083902B	8	11	39	135	0.035"	7
BXA085902B	8	11	59	135	0.035"	7
BXA087902B	8	11	79	135	0.035"	8
BXA092902B	9	13	29	135	0.035"	8
BXA093902B	9	13	39	135	0.035"	8
BXA095902B	9	13	59	135	0.035"	8
BXA097902B	9	13	79	135	0.035"	8
BXA102902B	10	13	29	135	0.035"	8
BXA103902B	10	13	39	135	0.035"	8
BXA105902B	10	13	59	135	0.035"	8
BXA107902B	10	13	79	135	0.035"	8
BXA112902B	11	16	29	135	0.035"	8
BXA113902B	11	16	39	135	0.035"	8
BXA115902B	11	16	59	135	0.035"	8
BXA117902B	11	16	79	135	0.035"	8
BXAL082902B	8	16	29	135	0.035"	7
BXAL083902B	8	16	39	135	0.035"	7
BXAL085902B	8	16	59	135	0.035"	8
BXAL087902B	8	16	79	135	0.035"	8

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

Products listed may not be available in all markets.

GORE, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2018–2020 W. L. Gore & Associates, Inc. AX1673-EN2 JANUARY 2020



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific) 800 437 8181 (United States)
00800 6334 4673 (Europe) 928 779 2771 (United States)

goremedical.com