



# GORE® VIATORR® TIPS 覆膜支架系

器械规格

# GORE® VIATORR® TIPS 覆膜支架系统

专为门静脉高压及其并发症的初治和修复治疗设计

产品编号	标示内径 (mm)	覆膜段长度 (cm)	裸段长度 (cm)	最大导丝直径 (英寸)	止血导引鞘 (Fr)	最大扩张球囊直径 (mm)
PTB084275W	8	4	2	≤ 0.038	10	8
PTB085275W	8	5	2	≤ 0.038	10	8
PTB086275W	8	6	2	≤ 0.038	10	8
PTB087275W	8	7	2	≤ 0.038	10	8
PTB088275W	8	8	2	≤ 0.038	10	8
PTB104275W	10	4	2	≤ 0.038	10	10
PTB105275W	10	5	2	≤ 0.038	10	10
PTB106275W	10	6	2	≤ 0.038	10	10
PTB107275W	10	7	2	≤ 0.038	10	10
PTB108275W	10	8	2	≤ 0.038	10	10

适用范围: GORE® VIATORR® TIPS 覆膜支架适用于初治和修复治疗门静脉高血压及其并发症, 如静脉曲张出血、胃病、顽固性腹水和 / 或肝性胸水。

 Consult Instructions  
for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

请参阅 [eifu.goremedical.com](http://eifu.goremedical.com) 上的“使用说明”, 以获取该产品销售市场的所有适用的适应症、警告、注意事项和禁忌症的完整说明。仅 Only

所列产品可能并非在所有市场中都有销售。

GORE, VIATORR and designs are trademarks of W. L. Gore & Associates.  
© 2020 W. L. Gore & Associates, Inc. AZ0872-SZH1 APRIL 2020



**W. L. GORE & ASSOCIATES, INC.**  
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific)  
1800 680 424 (Australia / New Zealand)  
00800 6334 4673 (Europe)  
800 437 8181 (United States)  
928 779 2771 (United States)

[goremedical.com](http://goremedical.com)