

GORE® VIABAHN® VBX
Balloon Expandable Endoprosthesis

CONFIGURATIONS



Together, improving life



Configurations

Catalogue Number	Stent Labeled/ Nominal Diameter (mm)	Max. Post-Dilated Stent Diameter (mm) ¹	Crimped Stent Length (mm)	Catheter Length (cm)	Guidewire Diameter	Introducer Sheath Size (Fr)
BXA051502E	5	8	15	135	0.035"	7
BXA051902E	5	8	19	135	0.035"	7
BXA052902E	5	8	29	135	0.035"	7
BXA053902E	5	8	39	135	0.035"	7
BXA055902E	5	8	59	135	0.035"	7
BXA057902E	5	8	79	135	0.035"	7
BXA061502E	6	8	15	135	0.035"	7
BXA061902E	6	8	19	135	0.035"	7
BXA062902E	6	8	29	135	0.035"	7
BXA063902E	6	8	39	135	0.035"	7
BXA065902E	6	8	59	135	0.035"	7
BXA067902E	6	8	79	135	0.035"	7
BXA071502E	7	11	15	135	0.035"	7
BXA071902E	7	11	19	135	0.035"	7
BXA072902E	7	11	29	135	0.035"	7
BXA073902E	7	11	39	135	0.035"	7
BXA075902E	7	11	59	135	0.035"	7
BXA077902E	7	11	79	135	0.035"	7
BXA082902E	8	11	29	135	0.035"	7
BXA083902E	8	11	39	135	0.035"	7
BXA085902E	8	11	59	135	0.035"	7
BXA087902E	8	11	79	135	0.035"	8
BXA092902E	9	13	29	135	0.035"	8
BXA093902E	9	13	39	135	0.035"	8
BXA095902E	9	13	59	135	0.035"	8
BXA097902E	9	13	79	135	0.035"	8
BXA102902E	10	13	29	135	0.035"	8
BXA103902E	10	13	39	135	0.035"	8
BXA105902E	10	13	59	135	0.035"	8
BXA107902E	10	13	79	135	0.035"	8
BXA112902E	11	16	29	135	0.035"	8
BXA113902E	11	16	39	135	0.035"	8
BXA115902E	11	16	59	135	0.035"	8
BXA117902E	11	16	79	135	0.035"	8
BXAL082902E	8	16	29	135	0.035"	7
BXAL083902E	8	16	39	135	0.035"	7
BXAL085902E	8	16	59	135	0.035"	8
BXAL087902E	8	16	79	135	0.035"	8

¹ Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. R_{only}

Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2021 W. L. Gore & Associates GmbH 21264550-EN AUGUST 2021

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86003 800 437 8181 928 779 2771

