

Gore Vascular Graft

# CONFIGURATIONS



*Together, improving life*





## GORE® PROPATEN® Vascular Graft

- H** Heparin, standard-walled, GORE® PROPATEN® Vascular Graft
- HR** Heparin, fixed ringed, standard-walled, GORE® PROPATEN® Vascular Graft
- HT** Heparin, thin-walled, GORE® PROPATEN® Vascular Graft
- HPT** Heparin, pediatric, thin-walled, GORE® PROPATEN® Vascular Graft
- HAX** Heparin axillobifemoral, removable ringed, standard-walled, GORE® PROPATEN® Vascular Graft
- IRH** Integrated ringed, standard-walled, heparin, GORE® PROPATEN® Vascular Graft
- IRTH** Integrated ringed, thin-walled, heparin, GORE® PROPATEN® Vascular Graft

## GORE-TEX® Vascular Graft

### GORE-TEX® STRETCH Vascular Graft



### GORE® INTERING® Vascular Graft

- IRS** GORE® INTERING® Vascular Graft, standard-walled, stretch
- IRST** GORE® INTERING® Vascular Graft, thin-walled, stretch
- R** Fixed ringed, standard-walled
- RR** Removable ringed, standard-walled
- RRT** Removable ringed, thin-walled
- S** Stretch, standard-walled
- SA** Stretch aortic, standard-walled
- SAX** Axillobifemoral, removable ringed, standard-walled
- SB** Stretch, bifurcated standard-walled limbs
- SBT** Stretch, bifurcated, thin-walled limbs
- SRRT** Stretch, removable ringed, thin-walled
- ST** Stretch, thin-walled
- L** Line (orientation)

### GORE® ACUSEAL Vascular Graft

- ECH** Early cannulation heparin, GORE® ACUSEAL Vascular Graft

KEY: PRODUCT NUMBERING SYSTEM (EXAMPLE ONLY)



**GORE-TEX® STRETCH Vascular Grafts** have a unique orientation line which aids in determining proper tension. When the longitudinal dots become dashes equal in length to the transverse marks, the graft is under proper tension.



**GORE® PROPATEN® Vascular Grafts**, which are based on Stretch technology, must be moderately tensioned prior to cutting the graft to length. The blue orientation markers can be used as an aid to determine moderate tension.



**Non-stretch GORE-TEX® Vascular Grafts** have a light blue dashed line, without transverse marks, and do not need to be tensioned before implantation.

## GORE® PROPATEN® Vascular Graft

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
<b>STANDARD-WALLED</b>							
H460045*	Yes	Stretch	Standard-walled	None	4–6	N/A	45
H470045*	Yes	Stretch	Standard-walled	None	4–7	N/A	45
H470080*	Yes	Stretch	Standard-walled	None	4–7	N/A	80
H050060	Yes	Stretch	Standard-walled	None	5	N/A	60
H060010	Yes	Stretch	Standard-walled	None	6	N/A	10
H060040	Yes	Stretch	Standard-walled	None	6	N/A	40
H060060	Yes	Stretch	Standard-walled	None	6	N/A	60
H060080	Yes	Stretch	Standard-walled	None	6	N/A	80
H080040	Yes	Stretch	Standard-walled	None	8	N/A	40
H080060	Yes	Stretch	Standard-walled	None	8	N/A	60
H080080	Yes	Stretch	Standard-walled	None	8	N/A	80
<b>STANDARD-WALLED FIXED RINGED</b>							
HR060545	Yes	Stretch	Standard-walled	Fixed ring	6	5	45
<b>STANDARD-WALLED INTEGRATED RINGED</b>							
IRH463845*	Yes	Stretch	Standard-walled	Integrated ring	4–6	38	45
IRH473845*	Yes	Stretch	Standard-walled	Integrated ring	4–7	38	45
IRH060520	Yes	Stretch	Standard-walled	Integrated ring	6	5	20
IRH061020	Yes	Stretch	Standard-walled	Integrated ring	6	10	20
IRH064040	Yes	Stretch	Standard-walled	Integrated ring	6	40	40
IRH060545	Yes	Stretch	Standard-walled	Integrated ring	6	5	45
<b>STANDARD-WALLED AXILLOBIFEMORAL REMOVABLE RINGED</b>							
HAX01	Yes	Stretch	Standard-walled	Removable ring	8 × 8	Fully ringed	70 × 40

\* Tapered

HAX02	Yes	Stretch	Standard-walled	Removable ring	8 × 8	Fully ringed	90 × 40
-------	-----	---------	-----------------	----------------	-------	--------------	---------

#### THIN-WALLED

HT060040	Yes	Stretch	Thin-walled	None	6	N/A	40
HT060080	Yes	Stretch	Thin-walled	None	6	N/A	80
HT070040	Yes	Stretch	Thin-walled	None	7	N/A	40
HT070080	Yes	Stretch	Thin-walled	None	7	N/A	80
HT080040	Yes	Stretch	Thin-walled	None	8	N/A	40
HT080080	Yes	Stretch	Thin-walled	None	8	N/A	80

#### THIN-WALLED REMOVABLE RINGED

HT057090	Yes	Stretch	Thin-walled	Removable ring	5	70	90
HT064050	Yes	Stretch	Thin-walled	Removable ring	6	40	50
HT063070	Yes	Stretch	Thin-walled	Removable ring	6	30	70
HT066070	Yes	Stretch	Thin-walled	Removable ring	6	60	70
HT066080	Yes	Stretch	Thin-walled	Removable ring	6	60	80
HT067090	Yes	Stretch	Thin-walled	Removable ring	6	70	90
HT074050	Yes	Stretch	Thin-walled	Removable ring	7	40	50
HT076080	Yes	Stretch	Thin-walled	Removable ring	7	60	80
HT083040	Yes	Stretch	Thin-walled	Removable ring	8	30	40
HT084050	Yes	Stretch	Thin-walled	Removable ring	8	40	50
HT087080	Yes	Stretch	Thin-walled	Removable ring	8	70	80
HT087090	Yes	Stretch	Thin-walled	Removable ring	8	70	90

#### THIN-WALLED INTEGRATED RINGED

IRTH064040	Yes	Stretch	Thin-walled	Integrated ring	6	40	40
IRTH066060	Yes	Stretch	Thin-walled	Integrated ring	6	60	60
IRTH068080	Yes	Stretch	Thin-walled	Integrated ring	6	80	80

IRTH074040	Yes	Stretch	Thin-walled	Integrated ring	7	40	40
IRTH076060	Yes	Stretch	Thin-walled	Integrated ring	7	60	60
IRTH078080	Yes	Stretch	Thin-walled	Integrated ring	7	80	80
IRTH084040	Yes	Stretch	Thin-walled	Integrated ring	8	40	40
IRTH086060	Yes	Stretch	Thin-walled	Integrated ring	8	60	60
IRTH088080	Yes	Stretch	Thin-walled	Integrated ring	8	80	80

## GORE® INTERING® Vascular Graft

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
<b>STANDARD-WALLED STRETCH</b>							
IRS46045LE*	No	Stretch	Standard-walled	Intering	4–6	38	45
IRS47045LE*	No	Stretch	Standard-walled	Intering	4–7	38	45
IRS06010020LE	No	Stretch	Standard-walled	Intering	6	10	20
IRS06040LE	No	Stretch	Standard-walled	Intering	6	40	40
IRS06005045LE	No	Stretch	Standard-walled	Intering	6	5	45
<b>THIN-WALLED STRETCH</b>							
IRST06040040LE	No	Stretch	Thin-walled	Intering	6	40	40
IRST06060060LE	No	Stretch	Thin-walled	Intering	6	60	60
IRST06080080LE	No	Stretch	Thin-walled	Intering	6	80	80
IRST07040040LE	No	Stretch	Thin-walled	Intering	7	40	40
IRST07060060LE	No	Stretch	Thin-walled	Intering	7	60	60
IRST07080080LE	No	Stretch	Thin-walled	Intering	7	80	80
IRST08040040LE	No	Stretch	Thin-walled	Intering	8	40	40
IRST08060060LE	No	Stretch	Thin-walled	Intering	8	60	60
IRST08080080LE	No	Stretch	Thin-walled	Intering	8	80	80

\* Tapered

## GORE-TEX® STRETCH Vascular Graft

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
<b>STANDARD-WALLED</b>							
S46045E*	No	Stretch	Standard-walled	None	4–6	N/A	45
S47030E*	No	Stretch	Standard-walled	None	4–7	N/A	30
S47045E*	No	Stretch	Standard-walled	None	4–7	N/A	45
S47055E*	No	Stretch	Standard-walled	None	4–7	N/A	55
S47070E*	No	Stretch	Standard-walled	None	4–7	N/A	70
S0504E	No	Stretch	Standard-walled	None	5	N/A	40
S0508E	No	Stretch	Standard-walled	None	5	N/A	80
S0602E	No	Stretch	Standard-walled	None	6	N/A	20
S0603E	No	Stretch	Standard-walled	None	6	N/A	30
S0604E	No	Stretch	Standard-walled	None	6	N/A	40
S0605E	No	Stretch	Standard-walled	None	6	N/A	50
S0607E	No	Stretch	Standard-walled	None	6	N/A	70
S0608E	No	Stretch	Standard-walled	None	6	N/A	80
S0704E	No	Stretch	Standard-walled	None	7	N/A	40
S0708E	No	Stretch	Standard-walled	None	7	N/A	80
S0804E	No	Stretch	Standard-walled	None	8	N/A	40
S0808E	No	Stretch	Standard-walled	None	8	N/A	80
S1004E	No	Stretch	Standard-walled	None	10	N/A	40
<b>THIN-WALLED</b>							
ST0407E	No	Stretch	Thin-walled	None	4	N/A	70
ST0507E	No	Stretch	Thin-walled	None	5	N/A	70
ST0601E	No	Stretch	Thin-walled	None	6	N/A	10
ST0602E	No	Stretch	Thin-walled	None	6	N/A	20

\* Tapered

ST0604E	No	Stretch	Thin-walled	None	6	N/A	40
ST0606E	No	Stretch	Thin-walled	None	6	N/A	60
ST0608E	No	Stretch	Thin-walled	None	6	N/A	80
ST0701E	No	Stretch	Thin-walled	None	7	N/A	10
ST0704E	No	Stretch	Thin-walled	None	7	N/A	40
ST0708E	No	Stretch	Thin-walled	None	7	N/A	80
ST0801E	No	Stretch	Thin-walled	None	8	N/A	10
ST0804E	No	Stretch	Thin-walled	None	8	N/A	40
ST0808E	No	Stretch	Thin-walled	None	8	N/A	80
ST1004E	No	Stretch	Thin-walled	None	10	N/A	40
ST1008E	No	Stretch	Thin-walled	None	10	N/A	80

#### THIN-WALLED REMOVABLE RINGED

SRRT05030040LE	No	Stretch	Thin-walled	Removable ring	5	30	40
SRRT05060070LE	No	Stretch	Thin-walled	Removable ring	5	60	70
SRRT06030040LE	No	Stretch	Thin-walled	Removable ring	6	30	40
SRRT06040050LE	No	Stretch	Thin-walled	Removable ring	6	40	50
SRRT06060080LE	No	Stretch	Thin-walled	Removable ring	6	60	80
SRRT06070080LE	No	Stretch	Thin-walled	Removable ring	6	70	80
SRRT08030040LE	No	Stretch	Thin-walled	Removable ring	8	30	40
SRRT08040050LE	No	Stretch	Thin-walled	Removable ring	8	40	50
SRRT08070080LE	No	Stretch	Thin-walled	Removable ring	8	70	80
SRRT08070100LE	No	Stretch	Thin-walled	Removable ring	8	70	100

#### STANDARD-WALLED LARGE DIAMETER

SI204E	No	Stretch	Standard-walled	None	12	N/A	40
SA1404E	No	Stretch	Standard-walled	None	14	N/A	40
SA1604E	No	Stretch	Standard-walled	None	16	N/A	40



SA1804E	No	Stretch	Standard-walled	None	18	N/A	40
SA2004E	No	Stretch	Standard-walled	None	20	N/A	40
SA2204E	No	Stretch	Standard-walled	None	22	N/A	40
SA2404E	No	Stretch	Standard-walled	None	24	N/A	40

#### STANDARD-WALLED AXILLOBIFEMORAL REMOVABLE RINGED

SAX01DE	No	Stretch	Standard-walled	Removable ring	8 x 8	Fully ringed	70 x 40
SAX02DE	No	Stretch	Standard-walled	Removable ring	8 x 8	Fully ringed	90 x 40

#### STANDARD-WALLED BIFURCATED

SB2201DE	No	Stretch	Standard-walled	None	22 x 11	N/A	40
SB2401DE	No	Stretch	Standard-walled	None	24 x 12	N/A	40

#### THIN-WALLED BIFURCATED

SBT1202DE	No	Stretch	Thin-walled	None	12 x 6	N/A	50
SBT1402DE	No	Stretch	Thin-walled	None	14 x 7	N/A	50
SBT1602DE	No	Stretch	Thin-walled	None	16 x 8	N/A	50
SBT1802DE	No	Stretch	Thin-walled	None	18 x 9	N/A	50
SBT2002DE	No	Stretch	Thin-walled	None	20 x 10	N/A	50

## GORE-TEX® Vascular Graft

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
------------------	-----------------------	--------------------	----------------	---------------------------	---------------	------------------------------	----------------------

#### THIN-WALLED REMOVABLE RINGED

RRT10070080LE	No	Non-stretch	Thin-walled	Removable ring	10	70	80
---------------	----	-------------	-------------	----------------	----	----	----

#### STANDARD-WALLED FIXED RINGED

R12030030E	No	Non-stretch	Standard-walled	Fixed ring	12	30	30
R14030030E	No	Non-stretch	Standard-walled	Fixed ring	14	30	30
R16030030E	No	Non-stretch	Standard-walled	Fixed ring	16	30	30
R20020020E	No	Non-stretch	Standard-walled	Fixed ring	20	20	20

## GORE® ACUSEAL Vascular Graft

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
ECH470045*	Yes	Non-stretch	N/A	None	4-7	N/A	45
ECH060040	Yes	Non-stretch	N/A	None	6	N/A	40

## GORE® PROPATEN® Vascular Graft configured for Pediatric Shunt

Catalogue number	CBAS® Heparin Surface	Stretch technology	Wall thickness	Radial support technology	Diameter (mm)	Radial supported length (cm)	Standard length (cm)
------------------	--------------------------	--------------------	----------------	---------------------------	---------------	------------------------------	----------------------

### STANDARD-WALLED

HPT030005	Yes	Stretch	Thin-walled	None	3	N/A	5
HPT030010	Yes	Stretch	Thin-walled	None	3	N/A	10
HPT030015	Yes	Stretch	Thin-walled	None	3	N/A	15
HPT350005	Yes	Stretch	Thin-walled	None	3.5	N/A	5
HPT350010	Yes	Stretch	Thin-walled	None	3.5	N/A	10
HPT350015	Yes	Stretch	Thin-walled	None	3.5	N/A	15
HPT040005	Yes	Stretch	Thin-walled	None	4	N/A	5
HPT040010	Yes	Stretch	Thin-walled	None	4	N/A	10
HPT040015	Yes	Stretch	Thin-walled	None	4	N/A	15
HPT050005	Yes	Stretch	Thin-walled	None	5	N/A	5
HPT050010	Yes	Stretch	Thin-walled	None	5	N/A	10
HPT050015	Yes	Stretch	Thin-walled	None	5	N/A	15
HPT060015	Yes	Stretch	Thin-walled	None	6	N/A	15

\* Tapered



 Consult Instructions  
for Use  
[eifu.goremedical.com](https://eifu.goremedical.com)

Refer to *Instructions for Use* at [eifu.goremedical.com](https://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx-Only

Products listed may not be available in all markets.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

GORE, GORE-TEX, *Together, improving life*, GORE ACUSEAL, INTERING, PROPATEN, STRETCH and designs are trademarks of W. L. Gore & Associates.  
© 2024 W. L. Gore & Associates GmbH 241404617-EN FEBRUARY 2024

---

**W. L. Gore & Associates, Inc.**  
[goremedical.com](https://goremedical.com)

**Asia Pacific** +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673  
**United States** Flagstaff, AZ 86004 800 437 8181 928 779 2771

