



► CONFIGURATIONS



PERFORMANCE
through innovation

CATALOGUE NUMBER	DIAMETER (mm)	LENGTH (mm)	CATHETER LENGTH (cm)	RECOMMENDED VESSEL DIAMETER	SHEATH COMPATIBILITY (Fr)	GUIDEWIRE COMPATIBILITY
PHB050301	5	30	80	4.0–4.7	6	0.035"
PHB050302	5	30	120	4.0–4.7	6	0.035"
PHB050401	5	40	80	4.0–4.7	6	0.035"
PHB050402	5	40	120	4.0–4.7	6	0.035"
PHB050601	5	60	80	4.0–4.7	6	0.035"
PHB050602	5	60	120	4.0–4.7	6	0.035"
PHB050801	5	80	80	4.0–4.7	6	0.035"
PHB050802	5	80	120	4.0–4.7	6	0.035"
PHB051001	5	100	80	4.0–4.7	6	0.035"
PHB051002	5	100	120	4.0–4.7	6	0.035"
PHB060301	6	30	80	4.8–5.5	6	0.035"
PHB060302	6	30	120	4.8–5.5	6	0.035"
PHB060401	6	40	80	4.8–5.5	6	0.035"
PHB060402	6	40	120	4.8–5.5	6	0.035"
PHB060601	6	60	80	4.8–5.5	6	0.035"
PHB060602	6	60	120	4.8–5.5	6	0.035"
PHB060801	6	80	80	4.8–5.5	6	0.035"
PHB060802	6	80	120	4.8–5.5	6	0.035"
PHB061001	6	100	80	4.8–5.5	6	0.035"
PHB061002	6	100	120	4.8–5.5	6	0.035"

Sizing, availability and pricing varies by country.
Please check with your Gore representative for availability.

CATALOGUE NUMBER	DIAMETER (mm)	LENGTH (mm)	CATHETER LENGTH (cm)	RECOMMENDED VESSEL DIAMETER	SHEATH COMPATIBILITY (Fr)	GUIDEWIRE COMPATIBILITY
PHB070301	7	30	80	5.6–6.5	6	0.035"
PHB070302	7	30	120	5.6–6.5	6	0.035"
PHB070401	7	40	80	5.6–6.5	6	0.035"
PHB070402	7	40	120	5.6–6.5	6	0.035"
PHB070601	7	60	80	5.6–6.5	6	0.035"
PHB070602	7	60	120	5.6–6.5	6	0.035"
PHB070801	7	80	80	5.6–6.5	6	0.035"
PHB070802	7	80	120	5.6–6.5	6	0.035"
PHB071001	7	100	80	5.6–6.5	6	0.035"
PHB071002	7	100	120	5.6–6.5	6	0.035"
PHB080301	8	30	80	6.6–7.5	7	0.035"
PHB080302	8	30	120	6.6–7.5	7	0.035"
PHB080401	8	40	80	6.6–7.5	7	0.035"
PHB080402	8	40	120	6.6–7.5	7	0.035"
PHB080601	8	60	80	6.6–7.5	7	0.035"
PHB080602	8	60	120	6.6–7.5	7	0.035"
PHB080801	8	80	80	6.6–7.5	7	0.035"
PHB080802	8	80	120	6.6–7.5	7	0.035"
PHB081001	8	100	80	6.6–7.5	7	0.035"
PHB081002	8	100	120	6.6–7.5	7	0.035"

Sizing, availability and pricing varies by country.
Please check with your Gore representative for availability.



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)	800.437.8181 (United States)
00800.6334.4673 (Europe)	928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.

GORE®, PERFORMANCE THROUGH INNOVATION, TIGRIS®, and designs are trademarks of W. L. Gore & Associates.
© 2011, 2014, 2016 W. L. Gore & Associates, Inc. AQ0664-EU3 SEPTEMBER 2016