

**25 cm**  
LONGEST  
LENGTH



— CONFIGURATIONS  
With PROPATEN Bioactive Surface



**VIABAHN<sup>®</sup>**  
ENDOPROSTHESIS

PROPATEN  
BIOACTIVE SURFACE

PERFORMANCE  
through innovation

CATALOGUE NUMBER	ENDOPROSTHESIS LABELED DIAMETER <sup>1</sup> (mm)	ENDOPROSTHESIS LENGTH <sup>1</sup> (cm)	CATHETER LENGTH <sup>2</sup> (cm)	RECOMMENDED VESSEL DIAMETER <sup>3</sup> (mm)	DEVICE PROFILE (Fr)	RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP (mm)
PAH050202	5	2.5	120	4.0 – 4.7	7	5
PAH050502	5	5.0	120	4.0 – 4.7	7	5
PAH051002	5	10.0	120	4.0 – 4.7	7	5
PAH051502	5	15.0	120	4.0 – 4.7	7	5
PAH052502	5	25.0	120	4.0 – 4.7	7	5
PAH060201	6	2.5	75	4.8 – 5.5	7	6
PAH060202	6	2.5	120	4.8 – 5.5	7	6
PAH060501	6	5.0	75	4.8 – 5.5	7	6
PAH060502	6	5.0	120	4.8 – 5.5	7	6
PAH061001	6	10.0	75	4.8 – 5.5	7	6
PAH061002	6	10.0	120	4.8 – 5.5	7	6
PAH061501	6	15.0	75	4.8 – 5.5	7	6
PAH061502	6	15.0	120	4.8 – 5.5	7	6
PAH062502	6	25.0	120	4.8 – 5.5	7	6
PAH070201	7	2.5	75	5.6 – 6.5	8	7
PAH070202	7	2.5	120	5.6 – 6.5	8	7
PAH070501	7	5.0	75	5.6 – 6.5	8	7
PAH070502	7	5.0	120	5.6 – 6.5	8	7
PAH071001	7	10.0	75	5.6 – 6.5	8	7
PAH071002	7	10.0	120	5.6 – 6.5	8	7
PAH071501	7	15.0	75	5.6 – 6.5	8	7
PAH071502	7	15.0	120	5.6 – 6.5	8	7
PAH072502	7	25.0	120	5.6 – 6.5	8	7
PAH080201	8	2.5	75	6.6 – 7.5	8	8
PAH080202	8	2.5	120	6.6 – 7.5	8	8
PAH080501	8	5.0	75	6.6 – 7.5	8	8
PAH080502	8	5.0	120	6.6 – 7.5	8	8
PAH081001	8	10.0	75	6.6 – 7.5	8	8
PAH081002	8	10.0	120	6.6 – 7.5	8	8
PAH081501	8	15.0	75	6.6 – 7.5	8	8
PAH081502	8	15.0	120	6.6 – 7.5	8	8
PAH082502	8	25.0	120	6.6 – 7.5	8	8

<sup>1</sup> Labeled device diameters and lengths are nominal.

<sup>2</sup> Ensure the guidewire is the appropriate size (see *Instructions for Use*) and has a length at least twice that of the delivery catheter.

<sup>3</sup> Recommended endoprosthesis compression within the vessel is approximately 5 – 20%.



CATALOGUE NUMBER	ENDOPROSTHESIS LABELED DIAMETER <sup>1</sup> (mm)	ENDOPROSTHESIS LENGTH <sup>1</sup> (cm)	CATHETER LENGTH <sup>2</sup> (cm)	RECOMMENDED VESSEL DIAMETER <sup>3</sup> (mm)	DEVICE PROFILE (Fr)	RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP <sup>5</sup> (mm)
PAH090502	9	5.0	120	7.6 – 8.5	9	9
PAH091002	9	10.0	120	7.6 – 8.5	9	9
PAH091502	9	15.0	120	7.6 – 8.5	9	9
PAH100502	10	5.0	120	8.6 – 9.5	11 <sup>4</sup>	10
PAH101002	10	10.0	120	8.6 – 9.5	11 <sup>4</sup>	10
PAH101502	10	15.0	120	8.6 – 9.5	11 <sup>4</sup>	10
PAH110502	11	5.0	120	9.6 – 10.5	11	12
PAH111002	11	10.0	120	9.6 – 10.5	11	12
PAH130502	13	5.0	120	10.6 – 12.0	12	14
PAH131002	13	10.0	120	10.6 – 12.0	12	14



<sup>1</sup> Labeled device diameters and lengths are nominal.

<sup>2</sup> Ensure the guidewire is the appropriate size (see *Instructions for Use*) and has a length at least twice that of the delivery catheter.

<sup>3</sup> Recommended endoprosthesis compression within the vessel is approximately 5 – 20%.

<sup>4</sup> The 10 mm diameter device is compatible with the following 10 Fr introducer sheaths: Cordis AVANTI® Sheath Introducer, Boston Scientific SUPER SHEATH Introducer Sheath, B. Braun INTRADYN Tear-Away Introducer Sheath.

<sup>5</sup> For the 11 and 13 mm diameter devices, balloon inflation pressure should not exceed 8 atm.



**W. L. GORE & ASSOCIATES, INC.**

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)

00800.6334.4673 (Europe)

800.437.8181 (United States)

928.779.2771 (United States)

[goremedical.com](http://goremedical.com)

 Consult Instructions for Use

Products listed may not be available in all markets.

AVANTI® is a trademark of Cordis Corporation.

SUPER SHEATH is a trademark of Boston Scientific.

INTRADYN is a trademark of B. Braun Medical, Inc.

GORE®, PERFORMANCE THROUGH INNOVATION, PROPATEN, TIP to HUB logo, VIABAHN®,

and designs are trademarks of W. L. Gore & Associates.

© 2007 – 2012 W. L. Gore & Associates, Inc. AR0326-EN1 JUNE 2012