



GORE[®] VIABAHN[®] VBX Balloon Expandable Endoprosthesis

CONFIGURATIONS



| Catalogue number | Stent labeled / nominal diameter (mm) | Max. post-dilated stent diameter (mm)* | Crimped stent length (mm) | Catheter length (cm) | Guidewire diameter | Introducer sheath size (Fr) |
|------------------|---------------------------------------|--|---------------------------|----------------------|--------------------|-----------------------------|
| BXA051502H | 5 | 8 | 15 | 135 | 0.035" | 7 |
| BXA051902H | 5 | 8 | 19 | 135 | 0.035" | 7 |
| BXA052902H | 5 | 8 | 29 | 135 | 0.035" | 7 |
| BXA053902H | 5 | 8 | 39 | 135 | 0.035" | 7 |
| BXA055902H | 5 | 8 | 59 | 135 | 0.035" | 7 |
| BXA057902H | 5 | 8 | 79 | 135 | 0.035" | 7 |
| BXA061502H | 6 | 8 | 15 | 135 | 0.035" | 7 |
| BXA061902H | 6 | 8 | 19 | 135 | 0.035" | 7 |
| BXA062902H | 6 | 8 | 29 | 135 | 0.035" | 7 |
| BXA063902H | 6 | 8 | 39 | 135 | 0.035" | 7 |
| BXA065902H | 6 | 8 | 59 | 135 | 0.035" | 7 |
| BXA067902H | 6 | 8 | 79 | 135 | 0.035" | 7 |
| BXA071502H | 7 | 11 | 15 | 135 | 0.035" | 7 |
| BXA071902H | 7 | 11 | 19 | 135 | 0.035" | 7 |
| BXA072902H | 7 | 11 | 29 | 135 | 0.035" | 7 |
| BXA073902H | 7 | 11 | 39 | 135 | 0.035" | 7 |
| BXA075902H | 7 | 11 | 59 | 135 | 0.035" | 7 |
| BXA077902H | 7 | 11 | 79 | 135 | 0.035" | 7 |
| BXA082902H | 8 | 11 | 29 | 135 | 0.035" | 7 |
| BXA083902H | 8 | 11 | 39 | 135 | 0.035" | 7 |
| BXA085902H | 8 | 11 | 59 | 135 | 0.035" | 7 |
| BXA087902H | 8 | 11 | 79 | 135 | 0.035" | 8 |
| BXA092902H | 9 | 13 | 29 | 135 | 0.035" | 8 |
| BXA093902H | 9 | 13 | 39 | 135 | 0.035" | 8 |
| BXA095902H | 9 | 13 | 59 | 135 | 0.035" | 8 |
| BXA097902H | 9 | 13 | 79 | 135 | 0.035" | 8 |
| BXA102902H | 10 | 13 | 29 | 135 | 0.035" | 8 |
| BXA103902H | 10 | 13 | 39 | 135 | 0.035" | 8 |
| BXA105902H | 10 | 13 | 59 | 135 | 0.035" | 8 |
| BXA107902H | 10 | 13 | 79 | 135 | 0.035" | 8 |
| BXA112902H | 11 | 16 [†] | 29 | 135 | 0.035" | 8 |
| BXA113902H | 11 | 16 [†] | 39 | 135 | 0.035" | 8 |
| BXA115902H | 11 | 16 [†] | 59 | 135 | 0.035" | 8 |
| BXA117902H | 11 | 16 [†] | 79 | 135 | 0.035" | 8 |
| BXAL082902H | 8 | 16 [†] | 29 | 135 | 0.035" | 7 |
| BXAL083902H | 8 | 16 [†] | 39 | 135 | 0.035" | 7 |
| BXAL085902H | 8 | 16 [†] | 59 | 135 | 0.035" | 8 |
| BXAL087902H | 8 | 16 [†] | 79 | 135 | 0.035" | 8 |

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

† Technical limit of the device. Expansion beyond 13 mm is outside of the indication – see *Instructions For Use*.

Products listed may not be available in all markets.

GORE, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2020 W. L. Gore & Associates, Inc. AY1620-EN1 JANUARY 2020



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific) 800 437 8181 (United States)
00800 6334 4673 (Europe) 928 779 2771 (United States)

goremedical.com