










## REFERENCE GUIDE

Shown below are GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement configurations for robotic and endoscopic surgical staplers:




### INTUITIVE® SUREFORM Stapler 60 Reloads

Tissue type	Thickness range (mm)	Package color indicator	GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement catalogue number	Closed staple height (mm*)
Vascular/thin	1.0–2.0		Not available for this reload	1.0
Regular	1.5–2.4		BSGX160GBE	1.5
Thick	2.0–3.3		BSGX160GBE	2.0
Very thick	2.3–4.0		BSGX160BKE	2.3

### ETHICON ECHELON ENDOPATH® Staple Line Reinforcement Reloads

Tissue type	Thickness range (mm)	Package color indicator	GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement catalogue number	Closed staple height (mm†)
Vascular/thin	1.0–2.0		BSGEC45E – Straight and Flex BSGEC60E – Straight BSGEC60AE – Flex	1.0
Regular	1.5–2.4		BSGEC45E – Straight and Flex BSGEC60E – Straight BSGEC60AE – Flex	1.5
Regular/thick	1.8–3.0		BSGEC45E – Straight and Flex BSGEC60E – Straight BSGEC60AE – Flex	1.8
Thick	2.0–3.3		BSGEC45E – Straight and Flex BSGEC60E – Straight BSGEC60AE – Flex	2.0
Very thick	2.3–4.0		BSGEC60E – Straight BSGEC60AE – Flex	2.3

### MEDTRONIC ENDO GIA Reinforced Reload with TRI-STAPLE Technology

Tissue type	Thickness range (mm)	Package color indicator	GORE® SEAMGUARD® Bioabsorbable Staple Line Reinforcement catalogue number	Closed staple height (mm‡)
Vascular/medium	N/A		Not available for this reload	0.88–1.5
Medium/thick	1.2–1.95		BSGTRI45PE BSGTRI60PE	1.5–2.25
Extra-thick	1.95–2.7		BSGTRI45BE BSGTRI60BE	2.25–3.0

\* Per the Stapler Company's *Instructions For Use*. INTUITIVE and SUREFORM are trademarks of Intuitive Surgical, Inc.

† Per the Stapler Company's *Instructions For Use*. ETHICON, ECHELON and ENDOPATH are trademarks of Ethicon, Inc.

‡ Per the Stapler Company's *Instructions For Use*. MEDTRONIC, ENDO GIA and TRI-STAPLE are trademarks of Medtronic, Inc.

 Consult Instructions  
for Use  
**eifu.goremedical.com**

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

ETHICON, ECHELON and ENDOPATH are trademarks of Ethicon, Inc.  
INTUITIVE and SUREFORM are trademarks of Intuitive Surgical, Inc.  
MEDTRONIC, ENDO GIA and TRI-STAPLE are trademarks of Medtronic, Inc  
GORE, *Together, improving life*, SEAMGUARD and designs are trademarks of W. L. Gore & Associates.  
© 2024 W. L. Gore & Associates, Inc. 24GP2043-EN01 JULY 2024

---

**W. L. Gore & Associates, Inc.**  
goremedical.com

**Asia Pacific** +65 67332882 **Europe** 00800 6334 4673  
**United States** Flagstaff, AZ 86004 800 437 8181 928 779 2771

