

# SIZING CHART



**GORE® ENFORM**  
Preperitoneal Biomaterial



**GORE® ENFORM**  
Intraperitoneal Biomaterial

*Together, improving life*



## GORE® ENFORM

### Preperitoneal Biomaterial

Catalogue number	Size	Minimum recommended trocar or trocar incision size*
GBWR0616	6 cm × 16 cm	10 mm trocar
GBWR0816	8 cm × 16 cm	11 mm trocar
GBWR1010	10 cm × 10 cm	12 mm trocar
GBWR1016	10 cm × 16 cm	12 mm trocar
GBWR1620	16 cm × 20 cm	19 mm trocar incision
GBWR2020	20 cm × 20 cm	24 mm trocar incision
GBWR2025	20 cm × 25 cm	24 mm trocar incision
GBWR2030	20 cm × 30 cm	24 mm trocar incision
GBWR2040	20 cm × 40 cm	24 mm trocar incision
GBWR2540	25 cm × 40 cm	24 mm trocar incision
GBWR3030	30 cm × 30 cm	28 mm trocar incision
GBWR3040	30 cm × 40 cm	28 mm trocar incision

## GORE® ENFORM

### Intraperitoneal Biomaterial

Catalogue number	Size	Minimum recommended trocar or trocar incision size*
GBFR0816	8 cm × 16 cm	11 mm trocar
GBFR1016	10 cm × 16 cm	12 mm trocar
GBFR1620	16 cm × 20 cm	21 mm trocar incision
GBFR2025	20 cm × 25 cm	24 mm trocar incision
GBFR2540	25 cm × 40 cm	28 mm trocar incision

\* Insert the device through an open incision, trocar incision, or introduce through a trocar. Excessive force should be avoided during trocar introduction to minimize potential for device damage. If GORE® ENFORM Biomaterial is not easily introduced down the trocar or trocar incision, replace the trocar with a larger sized trocar or augment the trocar incision size and retry.

Removing the trocar seal can help facilitate deployment.  
Ease of introduction may vary depending on rolled device size and trocar variability.



**6 cm × 16 cm**

Preperitoneal: GBWR0616

**8 cm × 16 cm**

Preperitoneal: GBWR0816

Intraperitoneal (with film): GBFR0816

**10 cm × 10 cm**

Preperitoneal: GBWR1010

**10 cm × 16 cm**

Preperitoneal: GBWR1016

Intraperitoneal (with film): GBFR1016

**16 cm × 20 cm**

Preperitoneal: GBWR1620

Intraperitoneal (with film): GBFR1620

**20 cm × 20 cm**

Preperitoneal: GBWR2020



**20 cm × 25 cm**

Preperitoneal: GBWR2025

Intraperitoneal (with film): GBFR2025

**20 cm × 30 cm**

Preperitoneal: GBWR2030

**20 cm × 40 cm**

Preperitoneal: GBWR2040

**25 cm x 40 cm**


Preperitoneal: GBWR2540

Intraperitoneal (with film): GBFR2540

**30 cm × 30 cm**


Preperitoneal: GBWR3030

**30 cm x 40 cm**  
Preperitoneal: GBWR3040



**GORE® ENFORM**

Preperitoneal Biomaterial





 Consult Instructions  
for Use  
[eifu.goremedical.com](https://eifu.goremedical.com)

Refer to *Instructions for Use* at [eifu.goremedical.com](https://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, ENFORM and designs are trademarks of W. L. Gore & Associates.  
© 2019, 2023 W. L. Gore & Associates, Inc. 23835202-EN FEBRUARY 2023

---

**W. L. Gore & Associates, Inc.**  
[goremedical.com](https://goremedical.com)

**Asia Pacific** +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673  
**United States** Flagstaff, AZ 86003 800 437 8181 928 779 2771

