



LOWER PROFILE VBX STENT GRAFT

Same endoprosthesis and trusted performance,
now with greater versatility.

**A lower delivery profile enables
a 1 Fr reduction on most sizes.**

Stent labeled/ nominal diameter (mm)	Stent length (mm)	Legacy profile (Fr)	Current profile (Fr)
5	15/19/29/39/59/79	7	6
6	15/19/29/39/59/79	7	6
7	15/19/29/39/59/79	7	6
8	29/39/59	7	7
8	79	8	7
8L	29/39	7	7
8L	59/79	8	8
9	29/39/59/79	8	7
10	29/39/59/79	8	8
11	29/39/59/79	8	8

NOW
6Fr
compatible

**GREATER
VERSATILITY**

All configurations available in .035" over-the-wire 80 cm and 135 cm catheter lengths.

The **longest** balloon expandable stent graft is now 6 Fr compatible and still offers the **widest** range of diameter adjustability.¹⁻³

Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*
5	8
6	8
7	11
8	11
8L	16†
9	13
10	13
11	16†

UP TO
79 mm
stent length¹⁻³



All configurations available in .035" over-the-wire 80 cm and 135 cm catheter lengths.

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

† Technical limit of the device as determined by in-vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the VBX FLEX Clinical Study and is outside of the approved indication — see *Instructions for Use*.

References

1. GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis [Instructions for Use]. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. MD169334.
2. LIFESTREAM® Balloon Expandable Vascular Covered Stent [Instructions for Use]. Tempe, AZ: Bard Peripheral Vascular, Inc; 2019. BAW1345700 Rev. 5 06/19.
3. iCast covered stent system [Instructions for Use]. Merrimack, NH: Atrium Medical Corporation; 2023. AW009603-EN Rev 11.

 Consult Instructions for Use
eifu.goremedical.com

INDICATIONS FOR USE IN THE U.S.: The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis is indicated for the treatment of de novo or restenotic lesions found in iliac arteries with reference vessel diameters ranging from 5 mm–13 mm and lesion lengths up to 110 mm, including lesions at the aortic bifurcation. **CONTRAINDICATIONS:** Do not use the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis in patients with known hypersensitivity to heparin, including those patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II. Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the market where this product is available. **Rx Only**

Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VBX FLEX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2023 W. L. Gore & Associates, Inc. 23968200-EN NOVEMBER 2023

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

