


NOW with Radiopaque Markers for Enhanced Visibility



Summary of changes:

- Addition of four gold radiopaque markers bonded to the graft at each end of endoprosthesis
- Only low profile (0.018" delivery system), 5–8 mm diameter devices incorporate this change
- Catalogue numbers will change for the 0.018" delivery system devices (PAJ configurations)

What has NOT changed:

- Delivery system and profile unchanged
- 0.035" delivery system devices (PAH configurations) do NOT have radiopaque markers on the endoprosthesis

W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface is sold in some markets as the GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface.

Products listed may not be available in all markets.

GORE®, PROPATEN, VIABAHN®, and designs are trademarks of W. L. Gore & Associates.
© 2017 W. L. Gore & Associates, Inc. AW1772-EN1 NOVEMBER 2017



PROPATEN
BIOACTIVE SURFACE