



# DEVICE SIZING GUIDE

## Device dimensions and sizing requirements

### GORE® TAG® Thoracic Branch Endoprosthesis with 8 mm portal

Aortic component and aortic component delivery system

Catalogue number	Device diameter (mm)	Intended aortic diameters (mm)	Internal portal diameter (mm)	Proximal segment length (mm)	Proximal covered length (mm)	Partially uncovered stent length (mm)	Overall device length (cm)	GORE® DrySeal Flex Introducer Sheath size (Fr)
TAC082110A	21	16–19.5	8	20	17	3	10	20
TAC082610A	26	19.5–24	8	20	16	4	10	20
TAC082615A	26	19.5–24	8	20	16	4	15	20
TAC082810A	28	22–26	8	20	16	4	10	22
TAC082815A	28	22–26	8	20	16	4	15	22
TAC082820A	28	22–26	8	20	16	4	20	22
TAC083115A	31	24–29	8	20	16	4	15	22
TAC083120A	31	24–29	8	20	16	4	20	22
TAC083415A	34	27–32	8	20	15	5	15	24
TAC083420A	34	27–32	8	20	15	5	20	24
TAC083715A	37	29–34	8	25	20	5	15	24
TAC083720A	37	29–34	8	25	20	5	20	24
TAC084015A	40	31–37	8	25	19	6	15	26
TAC084020A	40	31–37	8	25	19	6	20	26
TAC084515A	45	34–42	8	25	18.5	6.5	15	26
TAC084520A	45	34–42	8	25	18.5	6.5	20	26

### GORE® TAG® Thoracic Branch Endoprosthesis with 12 mm portal

Aortic component and aortic component delivery system

Catalogue number	Device diameter (mm)	Intended aortic diameters (mm)	Internal portal diameter (mm)	Proximal segment length (mm)	Proximal covered length (mm)	Partially uncovered stent length (mm)	Overall device length (cm)	GORE® DrySeal Flex Introducer Sheath size (Fr)
TAC123115A	31	24–29	12	40	36	4	15	22
TAC123120A	31	24–29	12	40	36	4	20	22
TAC123415A	34	27–32	12	40	35	5	15	24
TAC123420A	34	27–32	12	40	35	5	20	24
TAC123715A	37	29–34	12	40	35	5	15	24
TAC123720A	37	29–34	12	40	35	5	20	24
TAC124015A	40	31–37	12	40	34	6	15	26
TAC124020A	40	31–37	12	40	34	6	20	26
TAC124515A	45	34–42	12	40	33.5	6.5	15	26
TAC124520A	45	34–42	12	40	33.5	6.5	20	26

## Device dimensions and sizing requirements

### GORE® TAG® Thoracic Branch Endoprosthesis

Side branch component

Catalogue number	Device diameter (mm)	Intended branch vessel diameter (mm)	Portal segment diameter (mm)	Overall device length (cm)	Required branch vessel length (cm)
TSB080806A	8	6–7.5	8	6	3
TSB081006A	10	7.5–9	8	6	3
TSB081206A	12	9–11	8	6	3
TSB081506A	15	11–13	8	6	3
TSB081706A	17	13–15	8	6	3
TSB121506A	15	11–13	12	6	2.5
TSB121706A	17	13–15	12	6	2.5
TSB122006A	20	15–18	12	6	2.5

### GORE® TAG® Thoracic Branch Endoprosthesis

Aortic Extender and Aortic Extender delivery system

Catalogue number	Device diameter (mm)	Intended aortic diameters (mm)	Overall device length (mm)	Partially uncovered stent length (mm)	GORE® DrySeal Flex Introducer Sheath size (Fr)
TE2136A	21	16–19.5	3.6	3	20
TE2638A	26	19.5–24	3.8	4	20
TE2840A	28	22–26	4	4	22
TE3140A	31	24–29	4	4	22
TE3442A	34	27–32	4.2	5	24
TE3742A	37	29–34	4.2	5	24
TE4043A	40	31–37	4.3	6	26
TE4546A	45	34–42	4.6	6.5	26

 Consult Instructions for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

GORE, *Together, improving life* and designs are trademarks of W. L. Gore & Associates.

© 2022 W. L. Gore & Associates, Inc. 22467651-EN MAY 2022

W. L. Gore & Associates, Inc.  
[goremedical.com](http://goremedical.com)

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673  
 United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

