

POST-OP


# Performance from start to finish.

*On to better  
outcomes.*



## Advancing care through access.

*A new standard of flexibility to treat more challenging anatomies.*



**Deliver with Ease:** Hydrophilic coating and enhanced flexibility provide exceptional access to challenging anatomies and branch vessels

**Minimize Blood Loss:** Exclusive DrySeal valve enables introduction of multiple devices with proven hemostasis control

**Care for More Patients:** Optimized profile and new configurations provide tailored delivery options for a broader range of patient anatomy

**Complete Confidence:** Engineered for use with our endovascular portfolio



DRYSEAL

FLEX  
INTRODUCER SHEATH

## GORE® DrySeal Flex Introducer Sheath

CATALOGUE NUMBER	SHEATH SIZE (Fr)	MINIMUM SHEATH ID (mm)	NOMINAL SHEATH OD (mm)	WORKING LENGTH (cm)
DSF1233	12	4.0	4.7	33
DSF1245	12	4.0	4.7	45
DSF1433	14	4.7	5.3	33
DSF1533	15	5.0	5.6	33
DSF1633	16	5.3	6.1	33
DSF1833	18	6.0	6.7	33
DSF2033	20	6.7	7.5	33
DSF2065	20	6.7	7.5	65
DSF2233	22	7.3	8.2	33
DSF2265	22	7.3	8.2	65
DSF2433	24	8.0	8.8	33
DSF2465	24	8.0	8.8	65
DSF2633	26	8.7	9.5	33
DSF2665	26	8.7	9.5	65

Products listed may not be available in all markets.

# 12 Fr x 45 cm

*NOW AVAILABLE for use with the  
GORE® EXCLUDER® Iliac Branch Endoprosthesis*



**W. L. GORE & ASSOCIATES, INC.**  
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)  
00800.6334.4673 (Europe)  
800.437.8181 (United States)  
928.779.2771 (United States)

[goremedical.com](http://goremedical.com)