



CONFIGURATIONS

| Catalogue number | Description | Treatment range |
|------------------|-------------------------------------|-----------------|
| ASD27A | 27 mm GORE® CARDIOFORM ASD Occluder | 8–15 mm |
| ASD32A | 32 mm GORE® CARDIOFORM ASD Occluder | 13–20 mm |
| ASD37A | 37 mm GORE® CARDIOFORM ASD Occluder | 18–25 mm |
| ASD44A | 44 mm GORE® CARDIOFORM ASD Occluder | 23–30 mm |
| ASD48A | 48 mm GORE® CARDIOFORM ASD Occluder | 28–35 mm |

All device configurations have an 80 cm working length and 10–14 Fr working catheter depending on occluder size.

Refer to *Instructions for Use* at goremedical.com for a complete description of all contraindications, warnings, precautions, and adverse events.

Products listed may not be available in all markets.

GORE®, CARDIOFORM, TOGETHER, IMPROVING LIFE, and designs are trademarks of W. L. Gore & Associates, Inc.
© 2019 W. L. Gore & Associates, Inc. JUNE 2019

W. L. Gore & Associates, Inc.
Flagstaff, AZ 86004
+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States) goremedical.com



AY0861-ENT