

**GORE® EXCLUDER®**  
Endoprosthesis Family

# CONFIGURATIONS



*Together, improving life*



## GORE® EXCLUDER® AAA Endoprosthesis with C3® Delivery System

### Trunk-Ipsilateral Leg Endoprosthesis

Catalogue number	Intended aortic inner diameter (mm)	Aortic endoprosthesis diameter (mm)	Intended iliac inner diameter (mm)	Iliac endoprosthesis diameter (mm)	Overall device length (cm)	Recommended sheath size (Fr)
RLT231212	19–21	23	10–11	12	12	16
RLT231214	19–21	23	10–11	12	14	16
RLT231216	19–21	23	10–11	12	16	16
RLT231218	19–21	23	10–11	12	18	16
RLT231412	19–21	23	12–13.5	14.5	12	16
RLT231414	19–21	23	12–13.5	14.5	14	16
RLT231416	19–21	23	12–13.5	14.5	16	16
RLT231418	19–21	23	12–13.5	14.5	18	16
RLT261212	22–23	26	10–11	12	12	16
RLT261214	22–23	26	10–11	12	14	16
RLT261216	22–23	26	10–11	12	16	16
RLT261218	22–23	26	10–11	12	18	16
RLT261412	22–23	26	12–13.5	14.5	12	16
RLT261414	22–23	26	12–13.5	14.5	14	16
RLT261416	22–23	26	12–13.5	14.5	16	16
RLT261418	22–23	26	12–13.5	14.5	18	16
RLT281212	24–26	28.5	10–11	12	12	18
RLT281214	24–26	28.5	10–11	12	14	18
RLT281216	24–26	28.5	10–11	12	16	18
RLT281218	24–26	28.5	10–11	12	18	18
RLT281412	24–26	28.5	12–13.5	14.5	12	18
RLT281414	24–26	28.5	12–13.5	14.5	14	18
RLT281416	24–26	28.5	12–13.5	14.5	16	18
RLT281418	24–26	28.5	12–13.5	14.5	18	18
RLT311413	27–29	31	12–13.5	14.5	13	18
RLT311415	27–29	31	12–13.5	14.5	15	18
RLT311417	27–29	31	12–13.5	14.5	17	18
RLT351414	30–32	35	12–13.5	14.5	14	18
RLT351416	30–32	35	12–13.5	14.5	16	18
RLT351418	30–32	35	12–13.5	14.5	18	18

### Contralateral Leg Endoprosthesis

Catalogue number	Intended iliac inner diameter (mm)	Iliac endoprosthesis diameter (mm)	Contralateral Leg length (cm)	Recommended sheath size (Fr)
PLC121000	10–11	12	10	12
PLC121200	10–11	12	12	12
PLC121400	10–11	12	14	12
PLC141000	12–13.5	14.5	10	12
PLC141200	12–13.5	14.5	12	12
PLC141400	12–13.5	14.5	14	12
PLC161000*	13.5–14.5	16	9.5	12
PLC161200*	13.5–14.5	16	11.5	12
PLC161400*	13.5–14.5	16	13.5	12
PLC181000*	14.5–16.5	18	9.5	12
PLC181200*	14.5–16.5	18	11.5	12
PLC181400*	14.5–16.5	18	13.5	12
PLC201000*	16.5–18.5	20	9.5	12
PLC201200*	16.5–18.5	20	11.5	12
PLC201400*	16.5–18.5	20	13.5	12
PLC231000*	18.5–21.5	23	10	14
PLC231200*	18.5–21.5	23	12	14
PLC231400*	18.5–21.5	23	14	14
PLC271000*	21.5–25.0	27	10	15
PLC271200*	21.5–25.0	27	12	15
PLC271400*	21.5–25.0	27	14	15

### Aortic Extender Endoprosthesis

Catalogue number	Intended aortic inner diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended sheath size (Fr)
PLA230300	19–21	23	3.3	16
PLA260300	22–23	26	3.3	16
PLA280300	24–26	28.5	3.3	16
PLA320400	27–29	32	4.5	17
PLA360400	30–32	36	4.5	18

### Iliac Extender Endoprosthesis

Catalogue number	Intended iliac inner diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended sheath size (Fr)
PLL161007	8–9	10	7	12
PLL161207	10–11	12	7	12
PLL161407	12–13.5	14.5	7	12

## GORE® EXCLUDER® Iliac Branch Endoprosthesis

### Iliac Branch Component

Catalogue number	Intended external iliac vessel diameter (mm)	Diameter iliac end (mm)	Iliac branch component length (cm)	Recommended sheath size (Fr)
CEB231010A	6.5–9	10	10	16
CEB231210A	10–11	12	10	16
CEB231410A	12–13.5	14.5	10	16

### Internal Iliac Component

Catalogue number	Intended internal iliac vessel diameter (mm)	Diameter iliac end (mm)	Internal iliac component length (cm)	Recommended sheath size (Fr x cm)
HGB161007A	6.5–9	10	7	12 x 45
HGB161207A	10–11	12	7	12 x 45
HGB161407A	12–13.5	14.5	7	12 x 45

\* PLEASE NOTE: All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.

## Aortic Accessories

### GORE® DRYSEAL Flex Introducer Sheath

Catalogue number*	Sheath size (Fr)	Minimum sheath ID (mm)	Nominal sheath OD (mm)	Working length (cm)
GDSF1233	12	4.0	4.7	33
GDSF1245	12	4.0	4.7	45
GDSF1433	14	4.7	5.3	33
GDSF1533	15	5.0	5.6	33
GDSF1633	16	5.3	6.1	33
GDSF1833	18	6.0	6.7	33


\* GDSF catalogue numbers are not available in all regions; use DSF catalogue numbers instead where appropriate.

### GORE® Molding & Occlusion Balloon


Catalogue number	Size (Fr)	Balloon diameter (mm)	Catheter length (cm)
MOB37	10	10–37	90

# GORE® EXCLUDER® AAA Endoprosthesis featuring C3® Delivery System

## Aortic Extender and Trunk-Ipsilateral Leg Endoprosthesis





## Contralateral Leg Endoprosthesis<sup>1,2</sup>



## GORE® EXCLUDER® Iliac Branch Endoprosthesis System

### Internal Iliac Component Iliac Branch Component

### Iliac Extender Endoprosthesis<sup>3</sup>



\* 12 mm only available in select trunk sizes.  
† 14.5 mm only available in select trunk sizes.

1. All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.
2. The 27 mm x 10 cm Contralateral Leg Endoprosthesis does not have a radiopaque marker 6 cm above the distal end.
3. Traditional sizes of Iliac Extender (10, 12, 14.5 mm) are not compatible with large diameter (16, 18, 20, 23, 27 mm) Contralateral Leg Endoprostheses.


**Trunk-Ipsilateral Leg Endoprosthesis**

Catalogue number*	Intended aortic vessel diameter (mm)	Aortic endoprosthesis diameter (mm)	Intended iliac vessel diameter (mm)	Iliac endoprosthesis diameter (mm)	Overall device length (cm)	Recommended introducer sheath (Fr)
CXT201212	16–18	20	10–11	12	12	15
CXT201412	16–18	20	12–13.5	14.5	12	15
CXT231412	19–21	23	12–13.5	14.5	12	15
CXT261412	22–23	26	12–13.5	14.5	12	16
CXT281412	24–26	28.5	12–13.5	14.5	12	16
CXT321414	27–29	32	12–13.5	14.5	14	18
CXT361414	30–32	36	12–13.5	14.5	14	18

**Aortic Extender Endoprosthesis**

Catalogue number	Intended aortic inner diameter (mm)	Endoprosthesis diameter (mm)	Endoprosthesis length (cm)	Recommended sheath size (Fr)
CXA200005	16–18	20	4.5	15
CXA230005	19–21	23	4.5	15
CXA260005	22–23	26	4.5	15
CXA280005	24–26	28.5	4.5	16
CXA320005	27–29	32	4.5	18
CXA360005	30–32	36	4.5	18

**Aortic Extender and Trunk-Ipsilateral Leg Endoprosthesis**



\* Currently available commercially.

† 12 mm only available in the 20 mm GORE® EXCLUDER® Conformable AAA Endoprosthesis.

 Consult Instructions for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx only

Products listed may not be available in all markets.

GORE, *Together, improving life*, ACTIVE CONTROL, C3, DRYSEAL, EXCLUDER and designs are trademarks of W. L. Gore & Associates.

© 2021, 2024 W. L. Gore & Associates, Inc. 24AR1205-EN02 DECEMBER 2024

W. L. Gore & Associates, Inc.  
[goremedical.com](http://goremedical.com)

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673  
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

