



## THE TRUE HYBRID SOLUTION THAT PROMOTES HEALING

GORE® SYNECOR Biomaterial features Gore core technology to deliver high strength and lasting reliability.


- High-risk ventral repairs: High BMIs, comorbidities and recurrent hernias.
- Versatile, easy-to-use hybrid tri-layer solution.
- Tissue-building Gore 3D PGA:TMC web scaffold promotes rapid vascularity within seven days<sup>1</sup> and tissue ingrowth within one month.\*
- Intraperitoneal and preperitoneal configurations.
- Laparoscopic and open ventral hernia repair including robotic procedures.
- For bridging, where there is a need for permanent strength.

---

### MORE FOR YOUR PATIENTS

---

#### More size options now available



GORE® SYNECOR Preperitoneal Biomaterial

\* Data on file 2015; W.L. Gore & Associates, Inc; Flagstaff, AZ.



# Configurations



## GORE® SYNECOR Intraperitoneal Biomaterial

Catalogue number	Size
GKFC12	12 cm diameter*
GKfV1015	10 cm × 15 cm†
GKfV1520	15 cm × 20 cm†
GKfR2025	20 cm × 25 cm
GKfR2030	20 cm × 30 cm

### Trocar compatibility

GKfR2030 is designed to fit through a 12 mm trocar incision.

## GORE® SYNECOR Preperitoneal Biomaterial

Catalogue number	MORE SIZE OPTIONS	Minimum recommended trocar or trocar incision size
GKWC09	9 cm diameter* <b>Now available</b>	8 mm trocar
GKwV1015	10 cm x 15 cm†	10 mm trocar
GKwR1215	12 cm x 15 cm	11 mm trocar
GKwV1520	15 cm x 20 cm†	12 mm trocar
GKwR2025	20 cm x 25 cm	15 mm trocar
GKwR2030	20 cm x 30 cm	15 mm trocar‡
GKwR3030	30 cm x 30 cm <b>Now available</b>	18 mm trocar incision
GKwR3040	30 cm x 40 cm <b>Now available</b>	18 mm trocar incision
GKwR3535	35 cm x 35 cm <b>Now available</b>	21 mm trocar incision
GKwR4040	40 cm x 40 cm <b>Now available</b>	24 mm trocar incision
GKwR4050	40 cm x 50 cm <b>Now available</b>	24 mm trocar incision

\* Circle.

† Oval.

‡ Wetting recommended.

1. Crawford N. Assessment of Vascularity via Micro CT in Various Patch Devices. Flagstaff, AZ: W. L. Gore & Associates, Inc.; 2016. [Final study report]. 2344TL.

 Consult Instructions for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. FX only

Products listed may not be available in all markets.

GORE, *Together, improving life*, SYNECOR and designs are trademarks of W. L. Gore & Associates.

© 2021, 2022 W. L. Gore & Associates, Inc. 22435535-EN MARCH 2022

**W. L. Gore & Associates, Inc.**  
Flagstaff, AZ 86004

**Asia Pacific** +65 67332882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673  
**United States** 800 437 8181 928 779 2771 [goremedical.com](http://goremedical.com)

