

SIZING CHART



GORE® ENFORM
Preperitoneal Biomaterial



GORE® ENFORM
Intraperitoneal Biomaterial

Together, improving life



GORE® ENFORM

Preperitoneal Biomaterial

Catalogue number	Size	Minimum recommended trocar or trocar incision size*
GBWR0616	6 cm × 16 cm	10 mm trocar
GBWR0816	8 cm × 16 cm	11 mm trocar
GBWR1010	10 cm × 10 cm	12 mm trocar
GBWR1016	10 cm × 16 cm	12 mm trocar
GBWR1620	16 cm × 20 cm	19 mm trocar incision
GBWR2020	20 cm × 20 cm	24 mm trocar incision
GBWR2025	20 cm × 25 cm	24 mm trocar incision
GBWR2030	20 cm × 30 cm	24 mm trocar incision
GBWR2040	20 cm × 40 cm	24 mm trocar incision
GBWR2540	25 cm × 40 cm	24 mm trocar incision
GBWR3030	30 cm × 30 cm	28 mm trocar incision
GBWR3040	30 cm × 40 cm	28 mm trocar incision

GORE® ENFORM

Intraperitoneal Biomaterial

Catalogue number	Size	Minimum recommended trocar or trocar incision size*
GBFR0816	8 cm × 16 cm	11 mm trocar
GBFR1016	10 cm × 16 cm	12 mm trocar
GBFR1620	16 cm × 20 cm	21 mm trocar incision
GBFR2025	20 cm × 25 cm	24 mm trocar incision
GBFR2540	25 cm × 40 cm	28 mm trocar incision

* Insert the device through an open incision, trocar incision, or introduce through a trocar. Excessive force should be avoided during trocar introduction to minimize potential for device damage. If GORE® ENFORM Biomaterial is not easily introduced down the trocar or trocar incision, replace the trocar with a larger sized trocar or augment the trocar incision size and retry.

Removing the trocar seal can help facilitate deployment.
Ease of introduction may vary depending on rolled device size and trocar variability.



6 cm × 16 cm

Preperitoneal: GBWR0616

8 cm × 16 cm

Preperitoneal: GBWR0816

Intraperitoneal (with film): GBFR0816

10 cm × 10 cm

Preperitoneal: GBWR1010

10 cm × 16 cm

Preperitoneal: GBWR1016

Intraperitoneal (with film): GBFR1016

16 cm × 20 cm

Preperitoneal: GBWR1620

Intraperitoneal (with film): GBFR1620

20 cm × 20 cm

Preperitoneal: GBWR2020

20 cm × 25 cm

Preperitoneal: GBWR2025

Intraperitoneal (with film): GBFR2025

20 cm × 30 cm

Preperitoneal: GBWR2030

20 cm × 40 cm

Preperitoneal: GBWR2040

25 cm x 40 cm


Preperitoneal: GBWR2540

Intraperitoneal (with film): GBFR2540

30 cm × 30 cm


Preperitoneal: GBWR3030

30 cm x 40 cm
Preperitoneal: GBWR3040



GORE® ENFORM

Preperitoneal Biomaterial



 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, ENFORM and designs are trademarks of W. L. Gore & Associates.
© 2019, 2023 W. L. Gore & Associates, Inc. 23835202-EN FEBRUARY 2023

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86003 800 437 8181 928 779 2771

